

# lamp\_quick\_start

October 25, 2024

## 1 Quick Start

In this vignette we will demonstrate how to use LAMP python package. The input data and reference files are located in <https://github.com/wanchanglin/lamp/tree/master/examples/data>.

### 1.1 Setup

To use LAMP, the first step is to import some python libraries including LAMP.

```
[1]: import sqlite3
import pandas as pd
from lamp import anno, stats, utils
```

### 1.2 Data Loading

LAMP supports text files separated by comma (,) or tab (\t). The Microsoft's XLSX is also supported, provided that data set is in the first sheet.

Here we use a small example data set with tsv format. Load it into python and check its format:

```
[2]: d_data = "../data/df_pos_2.tsv"
data = pd.read_table(d_data, header=0, sep="\t")
data
```

```
[2]:
```

|     | name     | namecustom    | mz         | mzmin      | mzmax      | rt         | \ |
|-----|----------|---------------|------------|------------|------------|------------|---|
| 0   | M151T34  | M150.8867T34  | 150.886715 | 150.886592 | 150.886863 | 34.152700  |   |
| 1   | M151T40  | M151.0402T40  | 151.040235 | 151.040092 | 151.040350 | 39.838172  |   |
| 2   | M152T40  | M152.0436T40  | 152.043607 | 152.043451 | 152.043737 | 40.303700  |   |
| 3   | M153T34  | M152.8838T34  | 152.883824 | 152.883678 | 152.883959 | 34.174647  |   |
| 4   | M153T36  | M153.0195T36  | 153.019474 | 153.019331 | 153.019633 | 35.785847  |   |
| ..  | ...      | ...           | ...        | ...        | ...        | ...        |   |
| 395 | M283T339 | M283.2646T339 | 283.264583 | 283.264341 | 283.264809 | 338.763489 |   |
| 396 | M284T60  | M284.1953T60  | 284.195294 | 284.194939 | 284.195536 | 59.593561  |   |
| 397 | M284T108 | M284.2235T108 | 284.223499 | 284.223156 | 284.223692 | 108.406389 |   |
| 398 | M284T339 | M284.268T339  | 284.267962 | 284.267634 | 284.268204 | 338.725056 |   |
| 399 | M285T34  | M284.775T34   | 284.775031 | 284.774635 | 284.775287 | 34.079641  |   |

  

|   | rtmin     | rtmax     | npeaks | .  | ... | X210         | X209         | \ |
|---|-----------|-----------|--------|----|-----|--------------|--------------|---|
| 0 | 33.637595 | 35.465548 | 97     | 97 | ... | 4.224942e+06 | 3.946599e+06 |   |

|     |            |            |     |     |     |              |              |
|-----|------------|------------|-----|-----|-----|--------------|--------------|
| 1   | 37.556072  | 40.532315  | 95  | 95  | ... | 1.419062e+06 | 1.251606e+06 |
| 2   | 38.092678  | 40.909428  | 81  | 81  | ... | 1.203919e+05 | 9.970442e+04 |
| 3   | 33.637595  | 35.465548  | 98  | 98  | ... | 5.592065e+06 | 5.761380e+06 |
| 4   | 34.130244  | 36.287354  | 98  | 98  | ... | 7.284938e+06 | 1.083289e+07 |
| ..  | ...        | ...        | ... | ... | ... | ...          | ...          |
| 395 | 338.398380 | 339.165948 | 94  | 94  | ... | 3.509767e+05 | 4.117633e+05 |
| 396 | 58.844217  | 60.107058  | 59  | 59  | ... | NaN          | NaN          |
| 397 | 107.880510 | 108.971046 | 72  | 72  | ... | 7.477652e+04 | 7.482219e+04 |
| 398 | 338.268300 | 339.370098 | 84  | 84  | ... | 3.697604e+04 | 5.398264e+04 |
| 399 | 33.667172  | 35.198181  | 97  | 97  | ... | 3.439330e+06 | 3.359842e+06 |

  

|     | X208         | X207         | X206         | X205         | X204 \       |
|-----|--------------|--------------|--------------|--------------|--------------|
| 0   | 3.668948e+06 | 3.754321e+06 | 3.853724e+06 | 3.787350e+06 | 3.584464e+06 |
| 1   | 1.214826e+06 | 8.143028e+05 | 5.331963e+05 | 1.930928e+06 | 1.479001e+06 |
| 2   | 9.384000e+04 | 4.186335e+04 | NaN          | 2.115447e+05 | 1.285713e+05 |
| 3   | 5.845419e+06 | 5.576013e+06 | 5.552878e+06 | 6.132789e+06 | 5.891378e+06 |
| 4   | 1.140072e+07 | 8.220552e+06 | 9.255154e+06 | 7.648211e+06 | 7.723814e+06 |
| ..  | ...          | ...          | ...          | ...          | ...          |
| 395 | 3.948000e+05 | 4.338804e+05 | 5.335221e+05 | 6.224684e+05 | 7.009340e+05 |
| 396 | NaN          | NaN          | NaN          | 2.558004e+04 | 4.020517e+04 |
| 397 | 3.399667e+04 | 7.233564e+04 | 1.043879e+05 | 2.506785e+04 | 2.753769e+04 |
| 398 | 5.340109e+04 | 6.557698e+04 | 7.656575e+04 | 1.040606e+05 | 1.063727e+05 |
| 399 | 3.375577e+06 | 3.789056e+06 | 3.478506e+06 | 3.391588e+06 | 5.067802e+06 |

  

|     | X203         | X202         | X201         |
|-----|--------------|--------------|--------------|
| 0   | 3.499711e+06 | 3.623205e+06 | 4.145770e+06 |
| 1   | 1.076354e+06 | 9.293218e+05 | 5.298062e+05 |
| 2   | 9.389346e+04 | 7.163655e+04 | 4.916483e+04 |
| 3   | 5.418082e+06 | 5.036840e+06 | 5.733794e+06 |
| 4   | 5.571163e+06 | 5.362560e+06 | 9.259675e+06 |
| ..  | ...          | ...          | ...          |
| 395 | 3.005173e+05 | 3.133173e+05 | 8.204783e+05 |
| 396 | NaN          | 3.162670e+04 | 5.446684e+04 |
| 397 | NaN          | NaN          | NaN          |
| 398 | NaN          | 3.059370e+04 | 1.358056e+05 |
| 399 | 3.497546e+06 | 3.316025e+06 | 3.906000e+06 |

[400 rows x 110 columns]

This data set includes peak list and intensity data matrix. LAMP requires peak list's name, m/z value and retention time. User needs to indicate the locations of feature name, m/z value, retention time and starting points of data matrix from data. Here they are 1, 3, 6 and 11, respectively.

Load input data with `xlsx` format for LAMP:

```
[3]: cols = [1, 3, 6, 11]
      # d_data = "./data/df_pos_2.tsv"
      # df = anno.read_peak(d_data, cols, sep='\t')
```

```
d_data = "../data/df_pos_2.xlsx" # use xlsx file
df = anno.read_peak(d_data, cols)
df
```

```
[3]:
```

|     | name     | mz         | rt         | QC9          | QC5 \        |
|-----|----------|------------|------------|--------------|--------------|
| 0   | M151T34  | 150.886715 | 34.152700  | 3.664879e+06 | 3.735147e+06 |
| 1   | M151T40  | 151.040235 | 39.838172  | 7.406381e+05 | 7.524075e+05 |
| 2   | M152T40  | 152.043607 | 40.303700  | 6.105241e+04 | 5.335546e+04 |
| 3   | M153T34  | 152.883824 | 34.174647  | 5.141479e+06 | 5.496344e+06 |
| 4   | M153T36  | 153.019474 | 35.785847  | 5.336758e+06 | 5.558265e+06 |
| ..  | ...      | ...        | ...        | ...          | ...          |
| 395 | M283T339 | 283.264583 | 338.763489 | 7.330602e+05 | 8.243956e+05 |
| 396 | M284T60  | 284.195294 | 59.593561  | 2.310932e+04 | NaN          |
| 397 | M284T108 | 284.223499 | 108.406389 | 3.748444e+04 | 2.993283e+04 |
| 398 | M284T339 | 284.267962 | 338.725056 | 1.161886e+05 | 1.476514e+05 |
| 399 | M285T34  | 284.775031 | 34.079641  | 4.063268e+06 | 3.807148e+06 |

  

|     | QC4          | QC3          | QC26         | QC25         | QC24 \       |
|-----|--------------|--------------|--------------|--------------|--------------|
| 0   | 5.190263e+06 | 2.742966e+06 | 3.824723e+06 | 3.722932e+06 | 3.804188e+06 |
| 1   | NaN          | 6.429245e+05 | 1.167016e+06 | 1.175981e+06 | 1.122533e+06 |
| 2   | NaN          | NaN          | 6.875157e+04 | 7.807399e+04 | 8.943068e+04 |
| 3   | 8.335846e+06 | 3.860588e+06 | 5.316874e+06 | 5.988232e+06 | 5.844917e+06 |
| 4   | 1.118557e+07 | 6.876715e+06 | 9.967314e+06 | 9.073822e+06 | 9.328573e+06 |
| ..  | ...          | ...          | ...          | ...          | ...          |
| 395 | NaN          | 1.159506e+06 | 4.294760e+05 | 4.641813e+05 | 4.570657e+05 |
| 396 | NaN          | NaN          | 1.759336e+04 | 2.645392e+04 | 2.727266e+04 |
| 397 | NaN          | NaN          | 3.175596e+04 | 3.879604e+04 | 4.299529e+04 |
| 398 | NaN          | NaN          | NaN          | 6.753490e+04 | 5.436219e+04 |
| 399 | 4.645099e+06 | 2.232221e+06 | 4.576754e+06 | 4.533339e+06 | 4.559356e+06 |

  

|     | ... | X210         | X209         | X208         | X207 \       |
|-----|-----|--------------|--------------|--------------|--------------|
| 0   | ... | 4.224942e+06 | 3.946599e+06 | 3.668948e+06 | 3.754321e+06 |
| 1   | ... | 1.419062e+06 | 1.251606e+06 | 1.214826e+06 | 8.143028e+05 |
| 2   | ... | 1.203919e+05 | 9.970442e+04 | 9.384000e+04 | 4.186335e+04 |
| 3   | ... | 5.592065e+06 | 5.761380e+06 | 5.845419e+06 | 5.576013e+06 |
| 4   | ... | 7.284938e+06 | 1.083289e+07 | 1.140072e+07 | 8.220552e+06 |
| ..  | ... | ...          | ...          | ...          | ...          |
| 395 | ... | 3.509767e+05 | 4.117633e+05 | 3.948000e+05 | 4.338804e+05 |
| 396 | ... | NaN          | NaN          | NaN          | NaN          |
| 397 | ... | 7.477652e+04 | 7.482219e+04 | 3.399667e+04 | 7.233564e+04 |
| 398 | ... | 3.697604e+04 | 5.398264e+04 | 5.340109e+04 | 6.557698e+04 |
| 399 | ... | 3.439330e+06 | 3.359842e+06 | 3.375577e+06 | 3.789056e+06 |

  

|   | X206         | X205         | X204         | X203         | X202 \       |
|---|--------------|--------------|--------------|--------------|--------------|
| 0 | 3.853724e+06 | 3.787350e+06 | 3.584464e+06 | 3.499711e+06 | 3.623205e+06 |
| 1 | 5.331963e+05 | 1.930928e+06 | 1.479001e+06 | 1.076354e+06 | 9.293218e+05 |
| 2 | NaN          | 2.115447e+05 | 1.285713e+05 | 9.389346e+04 | 7.163655e+04 |

```

3    5.552878e+06  6.132789e+06  5.891378e+06  5.418082e+06  5.036840e+06
4    9.255154e+06  7.648211e+06  7.723814e+06  5.571163e+06  5.362560e+06
..    ...
395  5.335221e+05  6.224684e+05  7.009340e+05  3.005173e+05  3.133173e+05
396          NaN  2.558004e+04  4.020517e+04          NaN  3.162670e+04
397  1.043879e+05  2.506785e+04  2.753769e+04          NaN          NaN
398  7.656575e+04  1.040606e+05  1.063727e+05          NaN  3.059370e+04
399  3.478506e+06  3.391588e+06  5.067802e+06  3.497546e+06  3.316025e+06

```

```

X201
0    4.145770e+06
1    5.298062e+05
2    4.916483e+04
3    5.733794e+06
4    9.259675e+06
..    ...
395  8.204783e+05
396  5.446684e+04
397          NaN
398  1.358056e+05
399  3.906000e+06

```

[400 rows x 103 columns]

Data frame `df` now includes only `name`, `mz`, `rt` and intensity data matrix.

### 1.3 Metabolite Annotation

To perform metabolite annotation, users should provide their own reference file. Otherwise, LAMP will use its default reference file for annotation.

```

[4]: ref_path = "" # if empty, use default reference file for matching
# load reference library
cal_mass = False
ref = anno.read_ref(ref_path, calc=cal_mass)
ref

```

```

[4]:      compound_id molecular_formula      compound_name \
0          1638      C10Cl100      Chlordecone
1          38485      C10H10Br202      Dibromothymoquinone
2          32427      C10H10BrN02      Brofloxine (USAN/INN)
3          39834      C10H10Cl2N20      Fenmetozole (USAN)
4          10156      C10H10Cl203  4-(2,4-Dichlorophenoxy)butyric acid
...      ...      ...      ...
31639      80256      H5010P3      PPPi
31640      37374      H6N09P3      (Diphosphono)Aminophosphonic Acid
31641      32626      H9N204P      Ammonium phosphate (NF)
31642          735      HN03      Nitrate

```

|       |       |      |               |
|-------|-------|------|---------------|
| 31643 | 40762 | HN03 | Peroxynitrite |
|-------|-------|------|---------------|

|       | exact_mass |
|-------|------------|
| 0     | 485.683441 |
| 1     | 319.904755 |
| 2     | 254.989491 |
| 3     | 244.017018 |
| 4     | 248.000700 |
| ...   | ...        |
| 31639 | 257.909557 |
| 31640 | 256.925542 |
| 31641 | 132.029994 |
| 31642 | 62.995643  |
| 31643 | 62.995643  |

[31644 rows x 4 columns]

The reference file must have two columns: `molecular_formula` and `compound_name` (or `name`). The `exact_mass` is optional. if absent, LAMP will calculate 'exact\_mass' based on the NIST Atomic Weights and Isotopic Compositions for All Elements. If your reference file has `exact_mass` and you want to calculate it using NIST database, set `calc` as `True`. The `exact_mass` is used to match against a range of `mz`, controlled by `ppm`, in data frame `df`.

Another reference file is HMDB database for urine:

```
[5]: ref_path = "./data/hmdb_urine_v4_0_20200910_v1.tsv"
cal_mass = True # there is no exact mass in reference file, so calculate them
ref = anno.read_ref(ref_path, calc=cal_mass)
ref
```

```
[5]:
```

|      | id          | molecular_formula | molecular_name \                 |
|------|-------------|-------------------|----------------------------------|
| 0    | HMDB0000001 | C7H11N3O2         | 1-Methylhistidine                |
| 1    | HMDB0000002 | C3H10N2           | 1,3-Diaminopropane               |
| 2    | HMDB0000005 | C4H6O3            | 2-Ketobutyric acid               |
| 3    | HMDB0000008 | C4H8O3            | 2-Hydroxybutyric acid            |
| 4    | HMDB0000010 | C19H24O3          | 2-Methoxyestrone                 |
| ...  | ...         | ...               | ...                              |
| 1606 | HMDB0012308 | C8H8O3            | Vanillin                         |
| 1607 | HMDB0012322 | C10H8O            | 2-Naphthol                       |
| 1608 | HMDB0012325 | C5H10O5           | Arabinofuranose                  |
| 1609 | HMDB0012451 | C20H28O3          | all-trans-5,6-Epoxyretinoic acid |
| 1610 | HMDB0012467 | C15H13O9S         | (-)-Epicatechin sulfate          |

|   | inchi \   |
|---|---|
| 0 | InChI=1S/C7H11N3O2/c1-10-3-5(9-4-10)2-6(8)7(11... |
| 1 | InChI=1S/C3H10N2/c4-2-1-3-5/h1-5H2                |
| 2 | InChI=1S/C4H6O3/c1-2-3(5)4(6)7/h2H2,1H3,(H,6,7)   |
| 3 | InChI=1S/C4H8O3/c1-2-3(5)4(6)7/h3,5H,2H2,1H3,(... |

```

4      InChI=1S/C19H24O3/c1-19-8-7-12-13(15(19)5-6-18...
...
1606   InChI=1S/C8H8O3/c1-11-8-4-6(5-9)2-3-7(8)10/h2-...
1607   InChI=1S/C10H8O/c11-10-6-5-8-3-1-2-4-9(8)7-10/...
1608   InChI=1S/C5H10O5/c6-1-2-3(7)4(8)5(9)10-2/h2-9H...
1609   InChI=1S/C20H28O3/c1-15(8-6-9-16(2)14-17(21)22...
1610   InChI=1S/C15H14O9S/c16-9-3-8-5-13(24-25(20,21)...
```

|      | inchi_key                   | exact_mass |
|------|-----------------------------|------------|
| 0    | BRMWTNUJHUMWMS-LURJTMIESA-N | 169.085127 |
| 1    | XFNJVJPLKCPIBV-UHFFFAOYSA-N | 74.084398  |
| 2    | TYEYBOSBBBHJIV-UHFFFAOYSA-N | 102.031694 |
| 3    | AFENDNXGAFYKQO-VKHYHEASA-N  | 104.047344 |
| 4    | WHEUWNKSCXYKBU-QPWUGHHJSA-N | 300.172545 |
| ...  | ...                         | ...        |
| 1606 | MWOOGOJBHIARFG-UHFFFAOYSA-N | 152.047344 |
| 1607 | JWAZRIHNYRIHIV-UHFFFAOYSA-N | 144.057515 |
| 1608 | HMFHBZSHGGEWLO-HWQSCIPKSA-N | 150.052823 |
| 1609 | KEEHJLBAOLGBJZ-WEDZBJJJSA-N | 316.203845 |
| 1610 | WTXWEAXATVSZQX-AFYWNPRSA-M  | 369.028028 |

[1611 rows x 6 columns]

Next we use HMDB reference file for compounds match. Here function argument `ppm` is used to control the m/z matching tolerance(range).

```

[6]: ppm = 5.0
      match = anno.comp_match_mass(df, ppm, ref)
      match
```

| [6]: | id        | mz         | molecular_formula | molecular_name \          |
|------|-----------|------------|-------------------|---------------------------|
| 0    | M154T37   | 154.062402 | C8H10O3           | Hydroxytyrosol            |
| 1    | M164T119  | 164.046774 | C9H8O3            | Phenylpyruvic acid        |
| 2    | M164T119  | 164.046774 | C9H8O3            | m-Coumaric acid           |
| 3    | M164T119  | 164.046774 | C9H8O3            | 4-Hydroxycinnamic acid    |
| 4    | M164T119  | 164.046774 | C9H8O3            | 2-Hydroxycinnamic acid    |
| 5    | M164T233  | 164.046832 | C9H8O3            | Phenylpyruvic acid        |
| 6    | M164T233  | 164.046832 | C9H8O3            | m-Coumaric acid           |
| 7    | M164T233  | 164.046832 | C9H8O3            | 4-Hydroxycinnamic acid    |
| 8    | M164T233  | 164.046832 | C9H8O3            | 2-Hydroxycinnamic acid    |
| 9    | M164T53   | 164.046825 | C9H8O3            | Phenylpyruvic acid        |
| 10   | M164T53   | 164.046825 | C9H8O3            | m-Coumaric acid           |
| 11   | M164T53   | 164.046825 | C9H8O3            | 4-Hydroxycinnamic acid    |
| 12   | M164T53   | 164.046825 | C9H8O3            | 2-Hydroxycinnamic acid    |
| 13   | M167T35   | 167.021095 | C7H5NO4           | Quinolinic acid           |
| 14   | M173T36_3 | 173.104423 | C8H15NO3          | Hexanoylglycine           |
| 15   | M174T35   | 174.088395 | C8H14O4           | Suberic acid              |
| 16   | M181T36   | 181.060407 | C6H7N5O2          | 8-Hydroxy-7-methylguanine |

|    |         |            |             |                  |
|----|---------|------------|-------------|------------------|
| 17 | M212T39 | 212.067866 | C10H12O5    | Vanillactic acid |
| 18 | M276T36 | 276.077397 | C10H16N2O5S | Biotin sulfone   |

```

                                inchi \
0  InChI=1S/C8H10O3/c9-4-3-6-1-2-7(10)8(11)5-6/h1...
1  InChI=1S/C9H8O3/c10-8(9(11)12)6-7-4-2-1-3-5-7/...
2  InChI=1S/C9H8O3/c10-8-3-1-2-7(6-8)4-5-9(11)12/...
3  InChI=1S/C9H8O3/c10-8-4-1-7(2-5-8)3-6-9(11)12/...
4  InChI=1S/C9H8O3/c10-8-4-2-1-3-7(8)5-6-9(11)12/...
5  InChI=1S/C9H8O3/c10-8(9(11)12)6-7-4-2-1-3-5-7/...
6  InChI=1S/C9H8O3/c10-8-3-1-2-7(6-8)4-5-9(11)12/...
7  InChI=1S/C9H8O3/c10-8-4-1-7(2-5-8)3-6-9(11)12/...
8  InChI=1S/C9H8O3/c10-8-4-2-1-3-7(8)5-6-9(11)12/...
9  InChI=1S/C9H8O3/c10-8(9(11)12)6-7-4-2-1-3-5-7/...
10 InChI=1S/C9H8O3/c10-8-3-1-2-7(6-8)4-5-9(11)12/...
11 InChI=1S/C9H8O3/c10-8-4-1-7(2-5-8)3-6-9(11)12/...
12 InChI=1S/C9H8O3/c10-8-4-2-1-3-7(8)5-6-9(11)12/...
13 InChI=1S/C7H5NO4/c9-6(10)4-2-1-3-8-5(4)7(11)12...
14 InChI=1S/C8H15NO3/c1-2-3-4-5-7(10)9-6-8(11)12/...
15 InChI=1S/C8H14O4/c9-7(10)5-3-1-2-4-6-8(11)12/h...
16 InChI=1S/C6H7N5O2/c1-11-2-3(9-6(11)13)8-5(7)10...
17 InChI=1S/C10H12O5/c1-15-9-5-6(2-3-7(9)11)4-8(1...
18 InChI=1S/C10H16N2O5S/c13-8(14)4-2-1-3-7-9-6(5-...

```

|    | inchi_key                   | exact_mass | ppm_error |
|----|-----------------------------|------------|-----------|
| 0  | JUUBCHWRXWPFFH-UHFFFAOYSA-N | 154.06     | -3.84     |
| 1  | BTNMPGBKDVTSJY-UHFFFAOYSA-N | 164.05     | -3.47     |
| 2  | KKSDGJDHHZEWEP-SNAWJCMRSA-N | 164.05     | -3.47     |
| 3  | NGSWKAQJJWESNS-ZZXKWWIFSA-N | 164.05     | -3.47     |
| 4  | PMOWTIHVNWZYFI-AATRIKPKSA-N | 164.05     | -3.47     |
| 5  | BTNMPGBKDVTSJY-UHFFFAOYSA-N | 164.05     | -3.12     |
| 6  | KKSDGJDHHZEWEP-SNAWJCMRSA-N | 164.05     | -3.12     |
| 7  | NGSWKAQJJWESNS-ZZXKWWIFSA-N | 164.05     | -3.12     |
| 8  | PMOWTIHVNWZYFI-AATRIKPKSA-N | 164.05     | -3.12     |
| 9  | BTNMPGBKDVTSJY-UHFFFAOYSA-N | 164.05     | -3.16     |
| 10 | KKSDGJDHHZEWEP-SNAWJCMRSA-N | 164.05     | -3.16     |
| 11 | NGSWKAQJJWESNS-ZZXKWWIFSA-N | 164.05     | -3.16     |
| 12 | PMOWTIHVNWZYFI-AATRIKPKSA-N | 164.05     | -3.16     |
| 13 | GJAWHXHKYXBSV-UHFFFAOYSA-N  | 167.02     | -4.57     |
| 14 | UPCKIPHSMXJJOX-UHFFFAOYSA-N | 173.11     | -4.45     |
| 15 | TYFQFVWCELRYAO-UHFFFAOYSA-N | 174.09     | -4.67     |
| 16 | VHPXSVXJBWZORQ-UHFFFAOYSA-N | 181.06     | 2.39      |
| 17 | SVYIZYRTOYHQRE-UHFFFAOYSA-N | 212.07     | -2.86     |
| 18 | QPFQYMONYBAUCY-ZKWXMUHSA-N  | 276.08     | -2.16     |

match gives the compound matching results. LAMP also provides a mass adjust option by adduct library. You can provide your own adducts library otherwise LAMP uses its default adducts library.

The adducts library's format looks like:

```
[7]: add_path = './data/adducts_short.tsv'
lib_df = pd.read_csv(add_path, sep="\t")
lib_df
```

```
[7]:
```

|    | label                             | exact_mass | charge | ion_mode |
|----|-----------------------------------|------------|--------|----------|
| 0  | [M+H] <sup>+</sup>                | 1.007276   | 1      | pos      |
| 1  | [M+NH <sub>4</sub> ] <sup>+</sup> | 18.033826  | 1      | pos      |
| 2  | [M+Na] <sup>+</sup>               | 22.989221  | 1      | pos      |
| 3  | [M+Mg] <sup>+</sup>               | 23.984493  | 1      | pos      |
| 4  | [M+K] <sup>+</sup>                | 38.963158  | 1      | pos      |
| 5  | [M+Fe] <sup>+</sup>               | 55.934388  | 1      | pos      |
| 6  | [M+Cu] <sup>+</sup>               | 62.929049  | 1      | pos      |
| 7  | [M+2H] <sup>+</sup>               | 2.015101   | 1      | pos      |
| 8  | [M+3H] <sup>+</sup>               | 3.022926   | 1      | pos      |
| 9  | [M-H] <sup>-</sup>                | -1.007276  | 1      | neg      |
| 10 | [M+35Cl] <sup>-</sup>             | 34.969401  | 1      | neg      |
| 11 | [M+Formate] <sup>-</sup>          | 44.998203  | 1      | neg      |
| 12 | [M+Acetate] <sup>-</sup>          | 59.013853  | 1      | neg      |

The adducts library must have columns of label, exact\_mass, charge and ion\_mode.

We use this adducts file to adjust mass:

```
[8]: ion_mode = "pos"
# if empty, use default adducts library
add_path = './data/adducts_short.tsv'
lib_add = anno.read_lib(add_path, ion_mode)
lib_add
```

```
[8]:
```

|   | label                             | exact_mass | charge |
|---|-----------------------------------|------------|--------|
| 0 | [M+H] <sup>+</sup>                | 1.007276   | 1      |
| 1 | [M+NH <sub>4</sub> ] <sup>+</sup> | 18.033826  | 1      |
| 2 | [M+Na] <sup>+</sup>               | 22.989221  | 1      |
| 3 | [M+Mg] <sup>+</sup>               | 23.984493  | 1      |
| 4 | [M+K] <sup>+</sup>                | 38.963158  | 1      |
| 5 | [M+Fe] <sup>+</sup>               | 55.934388  | 1      |
| 6 | [M+Cu] <sup>+</sup>               | 62.929049  | 1      |
| 7 | [M+2H] <sup>+</sup>               | 2.015101   | 1      |
| 8 | [M+3H] <sup>+</sup>               | 3.022926   | 1      |

Now use function comp\_match\_mass\_add to match compounds:

```
[9]: match_1 = anno.comp_match_mass_add(df, ppm, ref, lib_add)
match_1
```

```
[9]:
```

|   | id      | mz         | molecular_formula   | molecular_name |
|---|---------|------------|---|----------------|
| 0 | M152T40 | 152.043607 | C <sub>5</sub> H <sub>8</sub> N <sub>2</sub> O <sub>2</sub> | Dihydrothymine |



|    |           |            |            |                                |
|----|-----------|------------|------------|--------------------------------|
| 1  | M154T37   | 154.062402 | C8H8O3     | p-Hydroxyphenylacetic acid     |
| 2  | M154T37   | 154.062402 | C8H8O3     | 3-Hydroxyphenylacetic acid     |
| 3  | M154T37   | 154.062402 | C8H8O3     | ortho-Hydroxyphenylacetic acid |
| 4  | M154T37   | 154.062402 | C8H8O3     | Mandelic acid                  |
| 5  | M154T37   | 154.062402 | C8H8O3     | 3-Cresotinic acid              |
| 6  | M154T37   | 154.062402 | C8H8O3     | 4-Hydroxy-3-methylbenzoic acid |
| 7  | M154T37   | 154.062402 | C8H8O3     | Vanillin                       |
| 8  | M157T35   | 157.036819 | C4H10N2O2  | 2,4-Diaminobutyric acid        |
| 9  | M157T35   | 157.036819 | C4H10N2O2  | L-2,4-diaminobutyric acid      |
| 10 | M167T35   | 167.021095 | C5H8N2O2   | Dihydrothymine                 |
| 11 | M174T35   | 174.088395 | C9H13NO    | Phenylpropanolamine            |
| 12 | M174T35   | 174.088395 | C10H14O    | Thymol                         |
| 13 | M174T35   | 174.088395 | C10H14O    | (S)-Carvone                    |
| 14 | M174T35   | 174.088395 | C8H12O4    | 2-Octenedioic acid             |
| 15 | M174T35   | 174.088395 | C8H12O4    | cis-4-Octenedioic acid         |
| 16 | M181T36   | 181.060407 | C8H8N2O3   | Nicotinuric acid               |
| 17 | M184T38   | 184.097942 | C10H13N2   | Nicotine imine                 |
| 18 | M185T39_2 | 185.082034 | C5H15N04P  | Phosphorylcholine              |
| 19 | M186T36   | 186.045606 | C6H14N2O   | N-Acetylputrescine             |
| 20 | M187T38   | 187.097642 | C5H15N04P  | Phosphorylcholine              |
| 21 | M193T40   | 193.050761 | C5H14N4    | Agmatine                       |
| 22 | M200T36   | 200.061328 | C7H16N2O   | N-Acetylcadaverine             |
| 23 | M201T39_1 | 201.051849 | C10H10O3   | 4-Methoxycinnamic acid         |
| 24 | M203T36_1 | 203.002108 | C9H9NO     | Indole-3-carbinol              |
| 25 | M212T39   | 212.067866 | C8H15NO3   | Hexanoylglycine                |
| 26 | M212T39   | 212.067866 | C10H10O5   | Vanilpyruvic acid              |
| 27 | M217T37_1 | 217.018279 | C10H11NO   | Tryptophol                     |
| 28 | M221T37   | 221.012328 | C9H11NO2   | L-Phenylalanine                |
| 29 | M223T38   | 223.008162 | C4H10NO6P  | O-Phosphothreonine             |
| 30 | M223T40   | 223.096863 | C12H14O4   | Monoisobutyl phthalic acid     |
| 31 | M226T44   | 226.128007 | C8H18N4O2  | Asymmetric dimethylarginine    |
| 32 | M226T44   | 226.128007 | C8H18N4O2  | Symmetric dimethylarginine     |
| 33 | M227T36   | 227.066175 | C9H10N2O5  | 3-Nitrotyrosine                |
| 34 | M229T38   | 229.069418 | C4H10N3O5P | Phosphocreatine                |
| 35 | M233T38   | 233.043479 | C8H10N4O2  | Caffeine                       |
| 36 | M245T44   | 245.045772 | C7H15N3O3  | Homocitrulline                 |
| 37 | M245T37_2 | 245.093315 | C13H18O2   | Ibuprofen                      |
| 38 | M249T38   | 249.038309 | C8H10N4O3  | 1,3,7-Trimethyluric acid       |
| 39 | M261T43   | 260.972975 | C10H7NO4   | Xanthurenic acid               |
| 40 | M269T37_2 | 269.088048 | C10H12N4O5 | Inosine                        |
| 41 | M275T168  | 275.201932 | C18H24O2   | Estradiol                      |
| 42 | M275T168  | 275.201932 | C18H24O2   | 17a-Estradiol                  |
| 43 | M277T181  | 277.217564 | C18H28O2   | 19-Norandrosterone             |
| 44 | M277T181  | 277.217564 | C18H28O2   | 19-Noretiocholanolone          |
| 45 | M278T71   | 278.148195 | C11H20N2O6 | Saccharopine                   |
| 46 | M279T233  | 279.233232 | C18H30O2   | alpha-Linolenic acid           |
| 47 | M279T233  | 279.233232 | C18H28O2   | 19-Norandrosterone             |

|    |          |            |            |                       |
|----|----------|------------|------------|-----------------------|
| 48 | M279T233 | 279.233232 | C18H28O2   | 19-Noretiocholanolone |
| 49 | M281T287 | 281.248903 | C18H32O2   | Linoleic acid         |
| 50 | M281T287 | 281.248903 | C18H30O2   | alpha-Linolenic acid  |
| 51 | M282T61  | 282.070271 | C10H14N2O6 | Ribothymidine         |
| 52 | M282T61  | 282.070271 | C10H14N2O6 | 3-Methyluridine       |
| 53 | M283T37  | 283.103695 | C11H14N4O5 | 1-Methylinosine       |

inchi \

```

0  InChI=1S/C5H8N2O2/c1-3-2-6-5(9)7-4(3)8/h3H,2H2...
1  InChI=1S/C8H8O3/c9-7-3-1-6(2-4-7)5-8(10)11/h1-...
2  InChI=1S/C8H8O3/c9-7-3-1-2-6(4-7)5-8(10)11/h1-...
3  InChI=1S/C8H8O3/c9-7-4-2-1-3-6(7)5-8(10)11/h1-...
4  InChI=1S/C8H8O3/c9-7(8(10)11)6-4-2-1-3-5-6/h1-...
5  InChI=1S/C8H8O3/c1-5-3-2-4-6(7(5)9)8(10)11/h2-...
6  InChI=1S/C8H8O3/c1-5-4-6(8(10)11)2-3-7(5)9/h2-...
7  InChI=1S/C8H8O3/c1-11-8-4-6(5-9)2-3-7(8)10/h2-...
8  InChI=1S/C4H10N2O2/c5-2-1-3(6)4(7)8/h3H,1-2,5-...
9  InChI=1S/C4H10N2O2/c5-2-1-3(6)4(7)8/h3H,1-2,5-...
10 InChI=1S/C5H8N2O2/c1-3-2-6-5(9)7-4(3)8/h3H,2H2...
11 InChI=1S/C9H13NO/c1-7(10)9(11)8-5-3-2-4-6-8/h2...
12 InChI=1S/C10H14O/c1-7(2)9-5-4-8(3)6-10(9)11/h4...
13 InChI=1S/C10H14O/c1-7(2)9-5-4-8(3)10(11)6-9/h4...
14 InChI=1S/C8H12O4/c9-7(10)5-3-1-2-4-6-8(11)12/h...
15 InChI=1S/C8H12O4/c9-7(10)5-3-1-2-4-6-8(11)12/h...
16 InChI=1S/C8H8N2O3/c11-7(12)5-10-8(13)6-2-1-3-9...
17 InChI=1S/C10H13N2/c1-12-7-3-5-10(12)9-4-2-6-11...
18 InChI=1S/C5H14NO4P/c1-6(2,3)4-5-10-11(7,8)9/h4...
19 InChI=1S/C6H14N2O/c1-6(9)8-5-3-2-4-7/h2-5,7H2,...
20 InChI=1S/C5H14NO4P/c1-6(2,3)4-5-10-11(7,8)9/h4...
21 InChI=1S/C5H14N4/c6-3-1-2-4-9-5(7)8/h1-4,6H2,(...
22 InChI=1S/C7H16N2O/c1-7(10)9-6-4-2-3-5-8/h2-6,8...
23 InChI=1S/C10H10O3/c1-13-9-5-2-8(3-6-9)4-7-10(1...
24 InChI=1S/C9H9NO/c11-6-7-5-10-9-4-2-1-3-8(7)9/h...
25 InChI=1S/C8H15NO3/c1-2-3-4-5-7(10)9-6-8(11)12/...
26 InChI=1S/C10H10O5/c1-15-9-5-6(2-3-7(9)11)4-8(1...
27 InChI=1S/C10H11NO/c12-6-5-8-7-11-10-4-2-1-3-9(...
28 InChI=1S/C9H11NO2/c10-8(9(11)12)6-7-4-2-1-3-5-...
29 InChI=1S/C4H10NO6P/c1-2(3(5)4(6)7)11-12(8,9)10...
30 InChI=1S/C12H14O4/c1-8(2)7-16-12(15)10-6-4-3-5...
31 InChI=1S/C8H18N4O2/c1-12(2)8(10)11-5-3-4-6(9)7...
32 InChI=1S/C8H18N4O2/c1-10-8(11-2)12-5-3-4-6(9)7...
33 InChI=1S/C9H10N2O5/c10-6(9(13)14)3-5-1-2-8(12)...
34 InChI=1S/C4H10N3O5P/c1-7(2-3(8)9)4(5)6-13(10,1...
35 InChI=1S/C8H10N4O2/c1-10-4-9-6-5(10)7(13)12(3)...
36 InChI=1S/C7H15N3O3/c8-5(6(11)12)3-1-2-4-10-7(9...
37 InChI=1S/C13H18O2/c1-9(2)8-11-4-6-12(7-5-11)10...
38 InChI=1S/C8H10N4O3/c1-10-4-5(9-7(10)14)11(2)8(...

```

39 InChI=1S/C10H7NO4/c12-7-3-1-2-5-8(13)4-6(10(14...  
 40 InChI=1S/C10H12N4O5/c15-1-4-6(16)7(17)10(19-4)...  
 41 InChI=1S/C18H24O2/c1-18-9-8-14-13-5-3-12(19)10...  
 42 InChI=1S/C18H24O2/c1-18-9-8-14-13-5-3-12(19)10...  
 43 InChI=1S/C18H28O2/c1-18-9-8-14-13-5-3-12(19)10...  
 44 InChI=1S/C18H28O2/c1-18-9-8-14-13-5-3-12(19)10...  
 45 InChI=1S/C11H20N2O6/c12-7(10(16)17)3-1-2-6-13-...  
 46 InChI=1S/C18H30O2/c1-2-3-4-5-6-7-8-9-10-11-12-...  
 47 InChI=1S/C18H28O2/c1-18-9-8-14-13-5-3-12(19)10...  
 48 InChI=1S/C18H28O2/c1-18-9-8-14-13-5-3-12(19)10...  
 49 InChI=1S/C18H32O2/c1-2-3-4-5-6-7-8-9-10-11-12-...  
 50 InChI=1S/C18H30O2/c1-2-3-4-5-6-7-8-9-10-11-12-...  
 51 InChI=1S/C10H14N2O6/c1-4-2-12(10(17)11-8(4)16)...  
 52 InChI=1S/C10H14N2O6/c1-11-6(14)2-3-12(10(11)17)...  
 53 InChI=1S/C11H14N4O5/c1-14-3-13-9-6(10(14)19)12...

|    | inchi_key                    | exact_mass | adduct              | ppm_error |
|----|------------------------------|------------|---------------------|-----------|
| 0  | NBAKTGXDIBVZOO-VKHYHEASA-N   | 152.04     | [M+Mg] <sup>+</sup> | 3.52      |
| 1  | XQXPVVBIMBYFF-UHFFFAOYSA-N   | 154.06     | [M+2H] <sup>+</sup> | -0.28     |
| 2  | FVMDYYGIDFPZAX-UHFFFAOYSA-N  | 154.06     | [M+2H] <sup>+</sup> | -0.28     |
| 3  | CCVYRRGZDBSHFU-UHFFFAOYSA-N  | 154.06     | [M+2H] <sup>+</sup> | -0.28     |
| 4  | IWYDHOAUDWTVEP-ZETCQYMHSA-N  | 154.06     | [M+2H] <sup>+</sup> | -0.28     |
| 5  | WHSXTWFYRGOBGO-UHFFFAOYSA-N  | 154.06     | [M+2H] <sup>+</sup> | -0.28     |
| 6  | LTFHNKUKQYVHDX-UHFFFAOYSA-N  | 154.06     | [M+2H] <sup>+</sup> | -0.28     |
| 7  | MWOOGOJBHIARFG-UHFFFAOYSA-N  | 154.06     | [M+2H] <sup>+</sup> | -0.28     |
| 8  | OGNSCSPNOLGXSM-UHFFFAOYSA-N  | 157.04     | [M+K] <sup>+</sup>  | -3.61     |
| 9  | OGNSCSPNOLGXSM-VKHYHEASA-N   | 157.04     | [M+K] <sup>+</sup>  | -3.61     |
| 10 | NBAKTGXDIBVZOO-VKHYHEASA-N   | 167.02     | [M+K] <sup>+</sup>  | -3.83     |
| 11 | DLNKOYKMWOXYQA-VXNVDRBHSA-N  | 174.09     | [M+Na] <sup>+</sup> | -3.10     |
| 12 | MGSRCZKZVOBKFT-UHFFFAOYSA-N  | 174.09     | [M+Mg] <sup>+</sup> | -3.23     |
| 13 | ULDHMXUKGWMISQ-VIFPVBQESA-N  | 174.09     | [M+Mg] <sup>+</sup> | -3.23     |
| 14 | BNTPVRGYUHHJFHN-HWKANZROSA-N | 174.09     | [M+2H] <sup>+</sup> | -1.52     |
| 15 | LQVYKEXVMZXOAH-UPHRSURJSA-N  | 174.09     | [M+2H] <sup>+</sup> | -1.52     |
| 16 | ZBSGKPYXQINNGF-UHFFFAOYSA-N  | 181.06     | [M+H] <sup>+</sup>  | -1.99     |
| 17 | GTQXYYYOJZZJHL-UHFFFAOYSA-N  | 184.10     | [M+Na] <sup>+</sup> | 4.61      |
| 18 | YHHSONZF0IEMCP-UHFFFAOYSA-O  | 185.08     | [M+H] <sup>+</sup>  | 4.80      |
| 19 | KLZGKIDSEJWEDW-UHFFFAOYSA-N  | 186.05     | [M+Fe] <sup>+</sup> | 3.25      |
| 20 | YHHSONZF0IEMCP-UHFFFAOYSA-O  | 187.10     | [M+3H] <sup>+</sup> | 4.52      |
| 21 | QYPPJABKJHAVHS-UHFFFAOYSA-N  | 193.05     | [M+Cu] <sup>+</sup> | -0.69     |
| 22 | RMOIHAKNOFHOE-UHFFFAOYSA-N   | 200.06     | [M+Fe] <sup>+</sup> | 3.39      |
| 23 | AFDXODALSZRGH-QPJXVBHSA-N    | 201.05     | [M+Na] <sup>+</sup> | -1.82     |
| 24 | IVYPNXXAYMYVSP-UHFFFAOYSA-N  | 203.00     | [M+Fe] <sup>+</sup> | -3.42     |
| 25 | UPCKIPHSXMXJOX-UHFFFAOYSA-N  | 212.07     | [M+K] <sup>+</sup>  | -2.29     |
| 26 | YGQHQTMZRPHIBB-UHFFFAOYSA-N  | 212.07     | [M+2H] <sup>+</sup> | -0.28     |
| 27 | MBBOMCVGYCRMEA-UHFFFAOYSA-N  | 217.02     | [M+Fe] <sup>+</sup> | -0.79     |
| 28 | COLNVLDPVHKLRT-QMMMGPOBSA-N  | 221.01     | [M+Fe] <sup>+</sup> | -4.70     |
| 29 | USRGIUJOYXOXQJ-GBXIJSLDSA-N  | 223.01     | [M+Mg] <sup>+</sup> | -4.06     |

|    |                              |        |                                   |       |
|----|------------------------------|--------|-----------------------------------|-------|
| 30 | RZJSUWQGFCNFS-UHFFFAOYSA-N   | 223.10 | [M+H] <sup>+</sup>                | 1.69  |
| 31 | YDGMGEXADBOMJ-LURJTMIESA-N   | 226.13 | [M+Mg] <sup>+</sup>               | 2.38  |
| 32 | HVPFXCBJHIIJGS-LURJTMIESA-N  | 226.13 | [M+Mg] <sup>+</sup>               | 2.38  |
| 33 | FBTSQILOGYXGMD-LURJTMIESA-N  | 227.07 | [M+H] <sup>+</sup>                | -0.32 |
| 34 | DRBBFCLWYRJSJZ-UHFFFAOYSA-N  | 229.07 | [M+NH <sub>4</sub> ] <sup>+</sup> | -0.94 |
| 35 | RYYVLZVUVIJVGH-UHFFFAOYSA-N  | 233.04 | [M+K] <sup>+</sup>                | -0.23 |
| 36 | XIGSAGMEBXLVJJ-YFKPBYRVSA-N  | 245.05 | [M+Fe] <sup>+</sup>               | 0.17  |
| 37 | HEFNNWSXXWATRW-UHFFFAOYSA-N  | 245.09 | [M+K] <sup>+</sup>                | -2.13 |
| 38 | BYXCUFUMGEBZDDI-UHFFFAOYSA-N | 249.04 | [M+K] <sup>+</sup>                | -0.56 |
| 39 | FBZONXHGGPHHIY-UHFFFAOYSA-N  | 260.97 | [M+Fe] <sup>+</sup>               | 4.13  |
| 40 | UGQMRVRMYASKQ-KQYNXXCUSA-N   | 269.09 | [M+H] <sup>+</sup>                | 0.01  |
| 41 | VOXZDWNPVJITMN-ZBRFXRBCSA-N  | 275.20 | [M+3H] <sup>+</sup>               | 5.00  |
| 42 | VOXZDWNPVJITMN-SFFUCWETSA-N  | 275.20 | [M+3H] <sup>+</sup>               | 5.00  |
| 43 | UOUIARGWRPHDBX-CQZDKXCPSA-N  | 277.22 | [M+H] <sup>+</sup>                | 4.90  |
| 44 | UOUIARGWRPHDBX-DHMOVHTBWSA-N | 277.22 | [M+H] <sup>+</sup>                | 4.90  |
| 45 | ZDGJAHTZVHVLOT-YUMQZZPRSA-N  | 278.15 | [M+2H] <sup>+</sup>               | 3.44  |
| 46 | DTOSIQBPPRVQHS-PDBXOOCHSA-N  | 279.23 | [M+H] <sup>+</sup>                | 4.93  |
| 47 | UOUIARGWRPHDBX-CQZDKXCPSA-N  | 279.23 | [M+3H] <sup>+</sup>               | 4.93  |
| 48 | UOUIARGWRPHDBX-DHMOVHTBWSA-N | 279.23 | [M+3H] <sup>+</sup>               | 4.93  |
| 49 | OYHQOLUKZRVRURQ-HZJYTTRNSA-N | 281.25 | [M+H] <sup>+</sup>                | 4.97  |
| 50 | DTOSIQBPPRVQHS-PDBXOOCHSA-N  | 281.25 | [M+3H] <sup>+</sup>               | 4.97  |
| 51 | DWRXFEITVBNRMK-JXOAFFINSA-N  | 282.07 | [M+Mg] <sup>+</sup>               | 2.10  |
| 52 | UTQUILVPBZEHTK-UHFFFAOYSA-N  | 282.07 | [M+Mg] <sup>+</sup>               | 2.10  |
| 53 | WJNGQIYEQLPJMN-IOSLPCCCSA-N  | 283.10 | [M+H] <sup>+</sup>                | -0.01 |

## 1.4 Correlation Analysis

Next step is correlation analysis, based on intensity data matrix along all peaks. All results are filtered by the correlation coefficient, p-values and retention time difference. That is: keep correlation results in an retention time differences/window (such as 1 second) with correlation coefficient larger than a threshold (such as 0.5) and their correlation p-values less than a threshold (such as 0.05).

LAMP supports two correlation methods, `pearson` and `spearman`. Also parameter `positive` allows user to select only positive correlation results, otherwise positive and negative correlations will be used.

Two functions, `_tic` and `_toc`, record the correlation computation time in seconds.

```
[10]: thres_rt = 1.0
      thres_corr = 0.5
      thres_pval = 0.05
      method = "spearman" # "pearson"
      positive = True
```

```
[11]: utils._tic()
      corr = stats.comp_corr_rt(df, thres_rt, thres_corr, thres_pval, method,
                               positive)
```

```
utils._toc()
corr
```

Elapsed time: 0.8097867965698242 seconds.

```
[11]:      name_a      name_b  r_value      p_value  rt_diff
0      M151T34      M153T34      0.80  1.267076e-23      0.02
1      M151T34      M155T34      0.71  1.752854e-16      0.20
2      M151T34      M161T34      0.78  1.869949e-21      0.14
3      M151T34      M163T34      0.69  3.239594e-15      0.20
4      M151T34      M167T35      0.51  5.776482e-08      0.73
...
1783  M283T34_1  M283T34_2      0.62  4.214876e-12      0.29
1784  M283T34_1      M285T34      0.82  5.937139e-26      0.08
1785  M283T34_2      M285T34      0.66  7.898957e-14      0.37
1786      M283T60      M284T60      0.86  1.033010e-29      0.15
1787  M283T339      M284T339      0.91  4.031333e-39      0.04
```

[1788 rows x 5 columns]

`corr` gives results of correlation coefficient(`r_value`), correlation p-values(`p_value`) and retention time difference(`rt_diff`).

Based on the correlation analysis, we can extract the groups and their sizes by:

```
[12]: # get correlation group and size
corr_df = stats.corr_grp_size(corr)
corr_df
```

```
[12]:      name  cor_grp_size      cor_grp
0      M217T35          52  M218T35::M219T34::M219T35::M221T34::M223T34::M...
1      M215T35          52  M216T35::M217T35::M218T35::M219T34::M219T35::M...
2      M216T35          52  M217T35::M218T35::M219T34::M219T35::M221T34::M...
3      M219T35          52  M221T34::M223T34::M225T35::M226T35::M229T34::M...
4      M218T35          51  M219T34::M219T35::M221T34::M223T34::M225T35::M...
..      ...          ...      ...
335     M266T66           1      M265T66
336     M193T40           1      M211T41
337     M191T233           1     M259T233
338     M278T181           1     M277T181
339     M277T71           1     M278T71
```

[340 rows x 3 columns]

## 1.5 Summarize Results

The final step gets the summary table in different format and save for the further analysis.

```
[13]: # get summary of metabolite annotation
sr, mr = anno.comp_summ(df, match)
```

This function combines peak table with compound matching results and returns two results in different formats. `sr` is single row results for each peak id in peak table `df`:

```
[14]: sr
```

```
[14]:
```

|     | name     | mz         | rt         | exact_mass | ppm_error | \ |
|-----|----------|------------|------------|------------|-----------|---|
| 0   | M151T34  | 150.886715 | 34.152700  | NaN        | NaN       |   |
| 1   | M151T40  | 151.040235 | 39.838172  | NaN        | NaN       |   |
| 2   | M152T40  | 152.043607 | 40.303700  | NaN        | NaN       |   |
| 3   | M153T34  | 152.883824 | 34.174647  | NaN        | NaN       |   |
| 4   | M153T36  | 153.019474 | 35.785847  | NaN        | NaN       |   |
| ..  | ...      | ...        | ...        | ...        | ...       |   |
| 395 | M283T61  | 283.068474 | 60.739869  | NaN        | NaN       |   |
| 396 | M284T108 | 284.223499 | 108.406389 | NaN        | NaN       |   |
| 397 | M284T339 | 284.267962 | 338.725056 | NaN        | NaN       |   |
| 398 | M284T60  | 284.195294 | 59.593561  | NaN        | NaN       |   |
| 399 | M285T34  | 284.775031 | 34.079641  | NaN        | NaN       |   |

|     | molecular_formula | molecular_name | inchi | inchi_key |
|-----|-------------------|----------------|-------|-----------|
| 0   | NaN               | NaN            | NaN   | NaN       |
| 1   | NaN               | NaN            | NaN   | NaN       |
| 2   | NaN               | NaN            | NaN   | NaN       |
| 3   | NaN               | NaN            | NaN   | NaN       |
| 4   | NaN               | NaN            | NaN   | NaN       |
| ..  | ...               | ...            | ...   | ...       |
| 395 | NaN               | NaN            | NaN   | NaN       |
| 396 | NaN               | NaN            | NaN   | NaN       |
| 397 | NaN               | NaN            | NaN   | NaN       |
| 398 | NaN               | NaN            | NaN   | NaN       |
| 399 | NaN               | NaN            | NaN   | NaN       |

[400 rows x 9 columns]

`mr` is multiple rows format if the match more than once from the reference file:

```
[15]: mr
```

```
[15]:
```

|     | name    | mz         | rt        | molecular_formula | molecular_name | inchi | \ |
|-----|---------|------------|-----------|-------------------|----------------|-------|---|
| 0   | M151T34 | 150.886715 | 34.152700 | NaN               | NaN            | NaN   |   |
| 1   | M151T40 | 151.040235 | 39.838172 | NaN               | NaN            | NaN   |   |
| 2   | M152T40 | 152.043607 | 40.303700 | NaN               | NaN            | NaN   |   |
| 3   | M153T34 | 152.883824 | 34.174647 | NaN               | NaN            | NaN   |   |
| 4   | M153T36 | 153.019474 | 35.785847 | NaN               | NaN            | NaN   |   |
| ..  | ...     | ...        | ...       | ...               | ...            | ...   |   |
| 404 | M283T61 | 283.068474 | 60.739869 | NaN               | NaN            | NaN   |   |

|     |          |            |            |     |     |     |
|-----|----------|------------|------------|-----|-----|-----|
| 405 | M284T108 | 284.223499 | 108.406389 | NaN | NaN | NaN |
| 406 | M284T339 | 284.267962 | 338.725056 | NaN | NaN | NaN |
| 407 | M284T60  | 284.195294 | 59.593561  | NaN | NaN | NaN |
| 408 | M285T34  | 284.775031 | 34.079641  | NaN | NaN | NaN |

|     | inchi_key | exact_mass | ppm_error |
|-----|-----------|------------|-----------|
| 0   | NaN       | NaN        | NaN       |
| 1   | NaN       | NaN        | NaN       |
| 2   | NaN       | NaN        | NaN       |
| 3   | NaN       | NaN        | NaN       |
| 4   | NaN       | NaN        | NaN       |
| ..  | ...       | ...        | ...       |
| 404 | NaN       | NaN        | NaN       |
| 405 | NaN       | NaN        | NaN       |
| 406 | NaN       | NaN        | NaN       |
| 407 | NaN       | NaN        | NaN       |
| 408 | NaN       | NaN        | NaN       |

[409 rows x 9 columns]

Now we merges single format results with correlation results:

```
[16]: # merge summery table with correlation analysis
res = anno.comp_summ_corr(sr, corr_df)
res
```

```
[16]:
```

|     | name     | mz         | rt         | exact_mass | ppm_error | \ |
|-----|----------|------------|------------|------------|-----------|---|
| 0   | M167T35  | 167.021095 | 34.882147  | 167.02     | -4.57     |   |
| 1   | M276T36  | 276.077397 | 36.385373  | 276.08     | -2.16     |   |
| 2   | M154T37  | 154.062402 | 37.183625  | 154.06     | -3.84     |   |
| 3   | M181T36  | 181.060407 | 35.734801  | 181.06     | 2.39      |   |
| 4   | M174T35  | 174.088395 | 35.001130  | 174.09     | -4.67     |   |
| ..  | ...      | ...        | ...        | ...        | ...       |   |
| 395 | M279T50  | 279.159930 | 50.055451  | NaN        | NaN       |   |
| 396 | M279T79  | 279.163910 | 78.758079  | NaN        | NaN       |   |
| 397 | M282T85  | 282.207859 | 84.719202  | NaN        | NaN       |   |
| 398 | M283T47  | 283.110871 | 46.822069  | NaN        | NaN       |   |
| 399 | M284T108 | 284.223499 | 108.406389 | NaN        | NaN       |   |

|     | molecular_formula | molecular_name            | \ |
|-----|-------------------|---------------------------|---|
| 0   | C7H5NO4           | Quinolinic acid           |   |
| 1   | C10H16N2O5S       | Biotin sulfone            |   |
| 2   | C8H10O3           | Hydroxytyrosol            |   |
| 3   | C6H7N5O2          | 8-Hydroxy-7-methylguanine |   |
| 4   | C8H14O4           | Suberic acid              |   |
| ..  | ...               | ...                       |   |
| 395 | NaN               | NaN                       |   |
| 396 | NaN               | NaN                       |   |

|     |     |     |
|-----|-----|-----|
| 397 | NaN | NaN |
| 398 | NaN | NaN |
| 399 | NaN | NaN |

  

|     | inchi   | \   |
|-----|---|-----|
| 0   | InChI=1S/C7H5NO4/c9-6(10)4-2-1-3-8-5(4)7(11)12... |     |
| 1   | InChI=1S/C10H16N2O5S/c13-8(14)4-2-1-3-7-9-6(5-... |     |
| 2   | InChI=1S/C8H10O3/c9-4-3-6-1-2-7(10)8(11)5-6/h1... |     |
| 3   | InChI=1S/C6H7N5O2/c1-11-2-3(9-6(11)13)8-5(7)10... |     |
| 4   | InChI=1S/C8H14O4/c9-7(10)5-3-1-2-4-6-8(11)12/h... |     |
| ..  | ...   |     |
| 395 |   | NaN |
| 396 |   | NaN |
| 397 |   | NaN |
| 398 |   | NaN |
| 399 |   | NaN |

  

|     | inchi_key                   | cor_grp_size | \ |
|-----|-----------------------------|--------------|---|
| 0   | GJAWHXHKYYXBSV-UHFFFAOYSA-N | 25.0         |   |
| 1   | QPFQYMONYBAUCY-ZKWXMUHSA-N  | 13.0         |   |
| 2   | JUUBCHWRXWPFFH-UHFFFAOYSA-N | 12.0         |   |
| 3   | VHPXSVXJBWZORQ-UHFFFAOYSA-N | 9.0          |   |
| 4   | TYFQFVWCELRYAO-UHFFFAOYSA-N | 9.0          |   |
| ..  | ...                         | ...          |   |
| 395 | NaN                         | NaN          |   |
| 396 | NaN                         | NaN          |   |
| 397 | NaN                         | NaN          |   |
| 398 | NaN                         | NaN          |   |
| 399 | NaN                         | NaN          |   |

  

|     | cor_grp  |
|-----|--|
| 0   | M171T34::M197T36::M209T34::M211T34::M213T34::M...  |
| 1   | M277T36_2::M278T36::M173T36_2::M186T36::M187T3...  |
| 2   | M155T38::M158T37_2::M164T36::M171T37_2::M173T3...  |
| 3   | M224T36::M225T35::M226T35::M227T36::M269T37_2::... |
| 4   | M211T34::M213T34::M219T34::M221T34::M229T35::M...  |
| ..  | ...  |
| 395 | NaN  |
| 396 | NaN  |
| 397 | NaN  |
| 398 | NaN  |
| 399 | NaN  |

[400 rows x 11 columns]

The result data frame **res** is re-arranged as four parts from top to bottom:

- 1st part: identified metabolites, satisfied with correlation analysis



- 2nd part: identified metabolites, not satisfied with correlation
- 3rd part: no identified metabolites, satisfied with correlation
- 4th part: no identified metabolites, not satisfied with correlation

The users should focus on the first part and perform their further analysis.

You can save all results in different forms, such as text format TSV or CSV. You can also save all results into a `sqlite3` database and use [DB Browser for SQLite](#) to view:

```
[17]: f_save = False          # here we do NOT save results
      db_out = "test.db"
      sr_out = "test_s.tsv"

[18]: if f_save:
      # save all results into a sqlite3 database
      conn = sqlite3.connect(db_out)
      df[["name", "mz", "rt"]].to_sql("peaklist",
                                     conn,
                                     if_exists="replace",
                                     index=False)

      corr_df.to_sql("corr_grp", conn, if_exists="replace", index=False)
      corr.to_sql("corr_pval_rt", conn, if_exists="replace", index=False)
      match.to_sql("match", conn, if_exists="replace", index=False)
      mr.to_sql("anno_mr", conn, if_exists="replace", index=False)
      res.to_sql("anno_sr", conn, if_exists="replace", index=False)

      conn.commit()
      conn.close()

      # save final results
      res.to_csv(sr_out, sep="\t", index=False)
```

## 1.6 End User Usages

For end users, LAMP provides two computation options: command line interface(CLI) and graphical user interface (GUI).

To use GUI, you need to open a terminal and type in:

```
$ lamp gui
```

To use CLI, open a terminal and type in command with required arguments, something like:

```
$ lamp cli \
  --sep "tab" \
  --input-data "./data/df_pos_3.tsv" \
  --col-idx "1, 2, 3, 4" \
  --add-path "" \
  --ref-path "" \
  --ion-mode "pos" \
  --cal-mass \
```

```

--thres-rt "1.0" \
--thres-corr "0.5" \
--thres-pval "0.05" \
--method "pearson" \
--positive \
--ppm "5.0" \
--save-db \
--save-mr \
--db-out "./res/test.db" \
--sr-out "./res/test_s.tsv" \
--mr-out "./res/test_m.tsv"

```

For the best practice, you can create a bash script `.sh` (Linux and MacOS) or Windows script `.bat` to contain these CLI arguments. Change parameters in these files each time when processing new data set.

For example, there are `lamp_cli.sh` and `lamp_cli.bat` in <https://github.com/wanchanglin/lamp/tree/master/examples>. You can run them and check the results in directory `examples/res`:

- For Linux and MacOS terminal:

```

$ chmod +x lamp_cli.sh
$ ./lamp_cli.sh

```

- For Windows terminal:

```

$ lamp_cli.bat

```