

lamp_quick_start

October 31, 2024

1 Quick Start

In this vignette we will demonstrate how to use LAMP python package. The input data and reference files are located in <https://github.com/wanchanglin/lamp/tree/master/examples/data>.

1.1 Setup

To use LAMP, the first step is to import some python libraries including LAMP.

```
[1]: import sqlite3
import pandas as pd
from lamp import anno, stats, utils
```

1.2 Data Loading

LAMP supports text files separated by comma (,) or tab (\t). The Microsoft's XLSX is also supported, provided that data set is in the first sheet.

Here we use a small example data set with tsv format. Load it into python and check its format:

```
[2]: d_data = "../data/df_pos_2.tsv"
data = pd.read_table(d_data, header=0, sep="\t")
data
```

```
[2]:
```

	name	namecustom	mz	mzmin	mzmax	rt	\
0	M151T34	M150.8867T34	150.886715	150.886592	150.886863	34.152700	
1	M151T40	M151.0402T40	151.040235	151.040092	151.040350	39.838172	
2	M152T40	M152.0436T40	152.043607	152.043451	152.043737	40.303700	
3	M153T34	M152.8838T34	152.883824	152.883678	152.883959	34.174647	
4	M153T36	M153.0195T36	153.019474	153.019331	153.019633	35.785847	
..	
395	M283T339	M283.2646T339	283.264583	283.264341	283.264809	338.763489	
396	M284T60	M284.1953T60	284.195294	284.194939	284.195536	59.593561	
397	M284T108	M284.2235T108	284.223499	284.223156	284.223692	108.406389	
398	M284T339	M284.268T339	284.267962	284.267634	284.268204	338.725056	
399	M285T34	M284.775T34	284.775031	284.774635	284.775287	34.079641	

	rtmin	rtmax	npeaks	X210	X209	\
0	33.637595	35.465548	97	97	...	4.224942e+06	3.946599e+06	

1	37.556072	40.532315	95	95	...	1.419062e+06	1.251606e+06
2	38.092678	40.909428	81	81	...	1.203919e+05	9.970442e+04
3	33.637595	35.465548	98	98	...	5.592065e+06	5.761380e+06
4	34.130244	36.287354	98	98	...	7.284938e+06	1.083289e+07
..
395	338.398380	339.165948	94	94	...	3.509767e+05	4.117633e+05
396	58.844217	60.107058	59	59	...	NaN	NaN
397	107.880510	108.971046	72	72	...	7.477652e+04	7.482219e+04
398	338.268300	339.370098	84	84	...	3.697604e+04	5.398264e+04
399	33.667172	35.198181	97	97	...	3.439330e+06	3.359842e+06

	X208	X207	X206	X205	X204 \
0	3.668948e+06	3.754321e+06	3.853724e+06	3.787350e+06	3.584464e+06
1	1.214826e+06	8.143028e+05	5.331963e+05	1.930928e+06	1.479001e+06
2	9.384000e+04	4.186335e+04	NaN	2.115447e+05	1.285713e+05
3	5.845419e+06	5.576013e+06	5.552878e+06	6.132789e+06	5.891378e+06
4	1.140072e+07	8.220552e+06	9.255154e+06	7.648211e+06	7.723814e+06
..
395	3.948000e+05	4.338804e+05	5.335221e+05	6.224684e+05	7.009340e+05
396	NaN	NaN	NaN	2.558004e+04	4.020517e+04
397	3.399667e+04	7.233564e+04	1.043879e+05	2.506785e+04	2.753769e+04
398	5.340109e+04	6.557698e+04	7.656575e+04	1.040606e+05	1.063727e+05
399	3.375577e+06	3.789056e+06	3.478506e+06	3.391588e+06	5.067802e+06

	X203	X202	X201
0	3.499711e+06	3.623205e+06	4.145770e+06
1	1.076354e+06	9.293218e+05	5.298062e+05
2	9.389346e+04	7.163655e+04	4.916483e+04
3	5.418082e+06	5.036840e+06	5.733794e+06
4	5.571163e+06	5.362560e+06	9.259675e+06
..
395	3.005173e+05	3.133173e+05	8.204783e+05
396	NaN	3.162670e+04	5.446684e+04
397	NaN	NaN	NaN
398	NaN	3.059370e+04	1.358056e+05
399	3.497546e+06	3.316025e+06	3.906000e+06

[400 rows x 110 columns]

This data set includes peak list and intensity data matrix. LAMP requires peak list's name, m/z value and retention time. User needs to indicate the locations of feature name, m/z value, retention time and starting points of data matrix from data. Here they are 1, 3, 6 and 11, respectively.

Load input data with `xlsx` format for LAMP:

```
[3]: ion_mode = "pos"
     cols = [1, 3, 6, 11]
     # d_data = "./data/df_pos_2.tsv"
```

```
# df = anno.read_peak(d_data, cols, sep='\t')
d_data = "./data/df_pos_2.xlsx" # use xlsx file
df = anno.read_peak(d_data, cols)
df
```

```
[3]:
```

	name	mz	rt	QC9	QC5 \
0	M151T34	150.886715	34.152700	3.664879e+06	3.735147e+06
1	M151T40	151.040235	39.838172	7.406381e+05	7.524075e+05
2	M152T40	152.043607	40.303700	6.105241e+04	5.335546e+04
3	M153T34	152.883824	34.174647	5.141479e+06	5.496344e+06
4	M153T36	153.019474	35.785847	5.336758e+06	5.558265e+06
..
395	M283T339	283.264583	338.763489	7.330602e+05	8.243956e+05
396	M284T60	284.195294	59.593561	2.310932e+04	NaN
397	M284T108	284.223499	108.406389	3.748444e+04	2.993283e+04
398	M284T339	284.267962	338.725056	1.161886e+05	1.476514e+05
399	M285T34	284.775031	34.079641	4.063268e+06	3.807148e+06

	QC4	QC3	QC26	QC25	QC24 \
0	5.190263e+06	2.742966e+06	3.824723e+06	3.722932e+06	3.804188e+06
1	NaN	6.429245e+05	1.167016e+06	1.175981e+06	1.122533e+06
2	NaN	NaN	6.875157e+04	7.807399e+04	8.943068e+04
3	8.335846e+06	3.860588e+06	5.316874e+06	5.988232e+06	5.844917e+06
4	1.118557e+07	6.876715e+06	9.967314e+06	9.073822e+06	9.328573e+06
..
395	NaN	1.159506e+06	4.294760e+05	4.641813e+05	4.570657e+05
396	NaN	NaN	1.759336e+04	2.645392e+04	2.727266e+04
397	NaN	NaN	3.175596e+04	3.879604e+04	4.299529e+04
398	NaN	NaN	NaN	6.753490e+04	5.436219e+04
399	4.645099e+06	2.232221e+06	4.576754e+06	4.533339e+06	4.559356e+06

	...	X210	X209	X208	X207 \
0	...	4.224942e+06	3.946599e+06	3.668948e+06	3.754321e+06
1	...	1.419062e+06	1.251606e+06	1.214826e+06	8.143028e+05
2	...	1.203919e+05	9.970442e+04	9.384000e+04	4.186335e+04
3	...	5.592065e+06	5.761380e+06	5.845419e+06	5.576013e+06
4	...	7.284938e+06	1.083289e+07	1.140072e+07	8.220552e+06
..
395	...	3.509767e+05	4.117633e+05	3.948000e+05	4.338804e+05
396	...	NaN	NaN	NaN	NaN
397	...	7.477652e+04	7.482219e+04	3.399667e+04	7.233564e+04
398	...	3.697604e+04	5.398264e+04	5.340109e+04	6.557698e+04
399	...	3.439330e+06	3.359842e+06	3.375577e+06	3.789056e+06

	X206	X205	X204	X203	X202 \
0	3.853724e+06	3.787350e+06	3.584464e+06	3.499711e+06	3.623205e+06
1	5.331963e+05	1.930928e+06	1.479001e+06	1.076354e+06	9.293218e+05

```

2          NaN  2.115447e+05  1.285713e+05  9.389346e+04  7.163655e+04
3  5.552878e+06  6.132789e+06  5.891378e+06  5.418082e+06  5.036840e+06
4  9.255154e+06  7.648211e+06  7.723814e+06  5.571163e+06  5.362560e+06
..          ...          ...          ...          ...
395  5.335221e+05  6.224684e+05  7.009340e+05  3.005173e+05  3.133173e+05
396          NaN  2.558004e+04  4.020517e+04          NaN  3.162670e+04
397  1.043879e+05  2.506785e+04  2.753769e+04          NaN          NaN
398  7.656575e+04  1.040606e+05  1.063727e+05          NaN  3.059370e+04
399  3.478506e+06  3.391588e+06  5.067802e+06  3.497546e+06  3.316025e+06

```

```

          X201
0  4.145770e+06
1  5.298062e+05
2  4.916483e+04
3  5.733794e+06
4  9.259675e+06
..          ...
395  8.204783e+05
396  5.446684e+04
397          NaN
398  1.358056e+05
399  3.906000e+06

```

[400 rows x 103 columns]

Data frame `df` now includes only `name`, `mz`, `rt` and intensity data matrix.

1.3 Metabolite Annotation

To perform metabolite annotation, users should provide their own reference file. Otherwise, LAMP will use its default reference file for annotation.

```

[4]: ref_path = "" # if empty, use default reference file for matching
      # load reference library
      cal_mass = False
      ref = anno.read_ref(ref_path, ion_mode=ion_mode, calc=cal_mass)
      ref

```

```

[4]:      compound_name  molecular_formula  monoisotopic_mass  \
34230      (-)-Salsoline      C11H15N02      193.110265
34231      (-)-trans-carveol      C10H160      152.120110
34232      (-)-ureidoglycolic acid      C3H6N2O4      134.032730
34233      (11R)-11-hydroperoxylinoleic acid      C18H32O4      312.230040
34234      (11Z,14Z)-eicosadienoylcarnitine      C27H49N04      451.366135
...          ...          ...          ...
83155  N(6),N(6),N(6)-trimethyl-L-lysine      C9H21N2O2+      189.160301
83156      nicotinic acid D-ribonucleotide      C11H15N09P+      336.048436
83157      phosphocholine      C5H15N04P+      184.073866

```

83158	S-adenosyl-L-methionine	C15H23N6O5S+	399.145060
83159	S-adenosylmethioninamine	C14H23N6O3S+	355.155232

	exact_mass	ion_type	ion_mode	\
34230	232.073425	[M+39K]+	positive	
34231	191.083270	[M+39K]+	positive	
34232	172.995890	[M+39K]+	positive	
34233	351.193200	[M+39K]+	positive	
34234	490.329295	[M+39K]+	positive	
...	
83155	189.159751	M+	positive	
83156	336.047886	M+	positive	
83157	184.073316	M+	positive	
83158	399.144510	M+	positive	
83159	355.154682	M+	positive	

	smiles	\
34230	<chem>CCc1cc2c(cc1O)CCN[C@H]2C</chem>	
34231	<chem>C=C(C)[C@@H]1CC=C(C)[C@@H](O)C1</chem>	
34232	<chem>NC(=O)N[C@@H](O)C(=O)O</chem>	
34233	<chem>CCCCC=CC(C=CCCCCCCC(=O)O)O</chem>	
34234	<chem>CCCCC/C=C\C/C=C\CCCCCCCCC(=O)OC(CC(=O)[O-])C[...]</chem>	
...	...	
83155	<chem>C[N+](C)(C)CCCC[C@H](N)C(=O)O</chem>	
83156	<chem>O=C(O)c1ccc[n+](c1)[C@@H]2O[C@H](COP(=O)(O)O)[C@@...</chem>	
83157	<chem>C[N+](C)(C)CCOP(=O)(O)O</chem>	
83158	<chem>C[S+](CC[C@H](N)C(=O)O)C[C@H]1O[C@@H](n2cnc3c(...</chem>	
83159	<chem>C[S+](CCCN)C[C@H]1O[C@@H](n2cnc3c(N)ncnc32)[C@...</chem>	

	inchikey	\
34230	YTPRLBGPZHUPO-ZETCQYMHSA-N	
34231	BAVONGHXFVOKBV-ZJUUDORDSA-N	
34232	NWZYCVIOKVITII-SFOWXEAESA-N	
34233	PLWDMWAXENHPLY-UHFFFAOYSA-N	
34234	OLZWDVKTGTGLC-UTJQPWESSA-N	
...	...	
83155	MXNRLFUSFKVQSK-QMMMGPOBSA-O	
83156	JOUIQRNQJGXQDC-ZYUZZMQFOSA-O	
83157	YHHSONZFOIEMCP-UHFFFAOYSA-O	
83158	MEFKEPWMEQBLKI-AIRLBKTGSA-O	
83159	ZUNBITIXDCPNDS-LSRJEVITSA-N	

	inchi	kegg_id	hmdb_id	\
34230	InChI=1S/C11H15N02/c1-7-9-6-11(14-2)10(13)5-8(...	C09640	-X-	
34231	InChI=1S/C10H16O/c1-7(2)9-5-4-8(3)10(11)6-9/h4...	C00964	-X-	
34232	InChI=1S/C3H6N2O4/c4-3(9)5-1(6)2(7)8/h1,6H,(H,...	C00603	HMDB0001005	
34233		-X-	-X-	-X-

```

34234 InChI=1S/C27H49N04/c1-5-6-7-8-9-10-11-12-13-14... -X- -X-
...
83155 InChI=1S/C9H20N2O2/c1-11(2,3)7-5-4-6-8(10)9(12... C03793 HMDB0001325
83156 InChI=1S/C11H14N09P/c13-8-7(5-20-22(17,18)19)2... C01185 -X-
83157 InChI=1S/C5H14N04P/c1-6(2,3)4-5-10-11(7,8)9/h4... C00588 HMDB0001565
83158 InChI=1S/C15H22N6O5S/c1-27(3-2-7(16)15(24)25)4... C00019 HMDB0001185
83159 InChI=1S/C14H23N6O3S/c1-24(4-2-3-15)5-8-10(21)... C01137 HMDB0000988

```

```

      chebi_id pubchem_id lipidmaps_id
34230 CHEBI:112      442356      -X-
34231 CHEBI:15389     -X-      -X-
34232 CHEBI:15412     439269     -X-
34233 CHEBI:134247    5230520     -X-
34234 CHEBI:73119     -X-      -X-
...
83155 CHEBI:17311     440120     -X-
83156 CHEBI:15763    53477721     -X-
83157 CHEBI:18132      1014     -X-
83158 CHEBI:15414    16757548     -X-
83159 CHEBI:15625     439415     -X-

```

[39150 rows x 14 columns]

The reference file must have two columns: `molecular_formula` and `compound_name` (or `name`). The `exact_mass` is optional. if absent, LAMP will calculate 'exact_mass' based on the NIST Atomic Weights and Isotopic Compositions for All Elements. If your reference file has `exact_mass` and you want to calculate it using NIST database, set `calc` as `True`. The `exact_mass` is used to match against a range of `mz`, controlled by `ppm`, in data frame `df`.

Another reference file is HMDB database for urine:

```

[5]: ref_path = "./data/hmdb_urine_v4_0_20200910_v1.tsv"
      cal_mass = True # there is no exact mass in reference file, so calculate them
      ref = anno.read_ref(ref_path, calc=cal_mass)
      ref

```

```

[5]:      id molecular_formula      compound_name \
0      HMDB0000001      C7H11N3O2      1-Methylhistidine
1      HMDB0000002      C3H10N2      1,3-Diaminopropane
2      HMDB0000005      C4H6O3      2-Ketobutyric acid
3      HMDB0000008      C4H8O3      2-Hydroxybutyric acid
4      HMDB0000010      C19H24O3      2-Methoxyestrone
...
1606   HMDB0012308      C8H8O3      Vanillin
1607   HMDB0012322      C10H8O      2-Naphthol
1608   HMDB0012325      C5H10O5      Arabinofuranose
1609   HMDB0012451      C20H28O3      all-trans-5,6-Epoxyretinoic acid
1610   HMDB0012467      C15H13O9S      (-)-Epicatechin sulfate

```

```

                                inchi \
0      InChI=1S/C7H11N3O2/c1-10-3-5(9-4-10)2-6(8)7(11...
1      InChI=1S/C3H10N2/c4-2-1-3-5/h1-5H2
2      InChI=1S/C4H6O3/c1-2-3(5)4(6)7/h2H2,1H3,(H,6,7)
3      InChI=1S/C4H8O3/c1-2-3(5)4(6)7/h3,5H,2H2,1H3,(...
4      InChI=1S/C19H24O3/c1-19-8-7-12-13(15(19)5-6-18...
...
1606   InChI=1S/C8H8O3/c1-11-8-4-6(5-9)2-3-7(8)10/h2-...
1607   InChI=1S/C10H8O/c11-10-6-5-8-3-1-2-4-9(8)7-10/...
1608   InChI=1S/C5H10O5/c6-1-2-3(7)4(8)5(9)10-2/h2-9H...
1609   InChI=1S/C20H28O3/c1-15(8-6-9-16(2)14-17(21)22...
1610   InChI=1S/C15H14O9S/c16-9-3-8-5-13(24-25(20,21)...
```

	inchi_key	exact_mass
0	BRMWTNUJHUMWMS-LURJTMIESA-N	169.085127
1	XFNJVJPLKCPIBV-UHFFFAOYSA-N	74.084398
2	TYEYBOSBBBHJIV-UHFFFAOYSA-N	102.031694
3	AFENDNXGAFYKQO-VKHYHEASA-N	104.047344
4	WHEUWNKSCXYKBU-QPWUGHHJSA-N	300.172545
...
1606	MWOOGOJBHIARFG-UHFFFAOYSA-N	152.047344
1607	JWAZRIHNYRIHIV-UHFFFAOYSA-N	144.057515
1608	HMFHBZSHGGEWLO-HWQSCIPKSA-N	150.052823
1609	KEEHJLBAOLGBJZ-WEDZBJJJSA-N	316.203845
1610	WTXWEAXATVSZQX-AFYWNPRSA-M	369.028028

[1611 rows x 6 columns]

Next we use HMDB reference file for compounds match. Here function argument ppm is used to control the m/z matching tolerance(range).

```

[6]: ppm = 5.0
      match = anno.comp_match_mass(df, ppm, ref)
      match
```

[6]:	id	mz	molecular_formula	compound_name
0	M154T37	154.062402	C8H10O3	Hydroxytyrosol
1	M164T119	164.046774	C9H8O3	Phenylpyruvic acid
2	M164T119	164.046774	C9H8O3	m-Coumaric acid
3	M164T119	164.046774	C9H8O3	4-Hydroxycinnamic acid
4	M164T119	164.046774	C9H8O3	2-Hydroxycinnamic acid
5	M164T233	164.046832	C9H8O3	Phenylpyruvic acid
6	M164T233	164.046832	C9H8O3	m-Coumaric acid
7	M164T233	164.046832	C9H8O3	4-Hydroxycinnamic acid
8	M164T233	164.046832	C9H8O3	2-Hydroxycinnamic acid
9	M164T53	164.046825	C9H8O3	Phenylpyruvic acid
10	M164T53	164.046825	C9H8O3	m-Coumaric acid

11	M164T53	164.046825	C9H8O3	4-Hydroxycinnamic acid
12	M164T53	164.046825	C9H8O3	2-Hydroxycinnamic acid
13	M167T35	167.021095	C7H5NO4	Quinolinic acid
14	M173T36_3	173.104423	C8H15NO3	Hexanoylglycine
15	M174T35	174.088395	C8H14O4	Suberic acid
16	M181T36	181.060407	C6H7N5O2	8-Hydroxy-7-methylguanine
17	M212T39	212.067866	C10H12O5	Vanillactic acid
18	M276T36	276.077397	C10H16N2O5S	Biotin sulfone

inchi \

```

0 InChI=1S/C8H10O3/c9-4-3-6-1-2-7(10)8(11)5-6/h1...
1 InChI=1S/C9H8O3/c10-8(9(11)12)6-7-4-2-1-3-5-7/...
2 InChI=1S/C9H8O3/c10-8-3-1-2-7(6-8)4-5-9(11)12/...
3 InChI=1S/C9H8O3/c10-8-4-1-7(2-5-8)3-6-9(11)12/...
4 InChI=1S/C9H8O3/c10-8-4-2-1-3-7(8)5-6-9(11)12/...
5 InChI=1S/C9H8O3/c10-8(9(11)12)6-7-4-2-1-3-5-7/...
6 InChI=1S/C9H8O3/c10-8-3-1-2-7(6-8)4-5-9(11)12/...
7 InChI=1S/C9H8O3/c10-8-4-1-7(2-5-8)3-6-9(11)12/...
8 InChI=1S/C9H8O3/c10-8-4-2-1-3-7(8)5-6-9(11)12/...
9 InChI=1S/C9H8O3/c10-8(9(11)12)6-7-4-2-1-3-5-7/...
10 InChI=1S/C9H8O3/c10-8-3-1-2-7(6-8)4-5-9(11)12/...
11 InChI=1S/C9H8O3/c10-8-4-1-7(2-5-8)3-6-9(11)12/...
12 InChI=1S/C9H8O3/c10-8-4-2-1-3-7(8)5-6-9(11)12/...
13 InChI=1S/C7H5NO4/c9-6(10)4-2-1-3-8-5(4)7(11)12...
14 InChI=1S/C8H15NO3/c1-2-3-4-5-7(10)9-6-8(11)12/...
15 InChI=1S/C8H14O4/c9-7(10)5-3-1-2-4-6-8(11)12/h...
16 InChI=1S/C6H7N5O2/c1-11-2-3(9-6(11)13)8-5(7)10...
17 InChI=1S/C10H12O5/c1-15-9-5-6(2-3-7(9)11)4-8(1...
18 InChI=1S/C10H16N2O5S/c13-8(14)4-2-1-3-7-9-6(5-...

```

	inchi_key	exact_mass	ppm_error
0	JUUBCHWRXWPFH-UHFFFAOYSA-N	154.06	-3.84
1	BTNMPGBKDVTSJY-UHFFFAOYSA-N	164.05	-3.47
2	KKSDGJDHHZEWEP-SNAWJCMRSA-N	164.05	-3.47
3	NGSWKAQJJWESNS-ZZXKWVIFSA-N	164.05	-3.47
4	PMOWTIHVNWZFYI-AATRIKPKSA-N	164.05	-3.47
5	BTNMPGBKDVTSJY-UHFFFAOYSA-N	164.05	-3.12
6	KKSDGJDHHZEWEP-SNAWJCMRSA-N	164.05	-3.12
7	NGSWKAQJJWESNS-ZZXKWVIFSA-N	164.05	-3.12
8	PMOWTIHVNWZFYI-AATRIKPKSA-N	164.05	-3.12
9	BTNMPGBKDVTSJY-UHFFFAOYSA-N	164.05	-3.16
10	KKSDGJDHHZEWEP-SNAWJCMRSA-N	164.05	-3.16
11	NGSWKAQJJWESNS-ZZXKWVIFSA-N	164.05	-3.16
12	PMOWTIHVNWZFYI-AATRIKPKSA-N	164.05	-3.16
13	GJAWHXHKYYXBSV-UHFFFAOYSA-N	167.02	-4.57
14	UPCKIPHSMXJOX-UHFFFAOYSA-N	173.11	-4.45
15	TYFQFVWCELRYAO-UHFFFAOYSA-N	174.09	-4.67

16	VHPXSVXJBWZORQ-UHFFFAOYSA-N	181.06	2.39
17	SVYIZYRTOYHQRE-UHFFFAOYSA-N	212.07	-2.86
18	QPFQYMONYBAUCY-ZKWXMUHSA-N	276.08	-2.16

`match` gives the compound matching results. LAMP also provides a mass adjust option by adduct library. You can provide your own adducts library otherwise LAMP uses its default adducts library.

The adducts library's format looks like:

```
[7]: add_path = './data/adducts_short.tsv'
lib_df = pd.read_csv(add_path, sep="\t")
lib_df
```

```
[7]:
```

	label	exact_mass	charge	ion_mode
0	[M+H] ⁺	1.007276	1	pos
1	[M+NH ₄] ⁺	18.033826	1	pos
2	[M+Na] ⁺	22.989221	1	pos
3	[M+Mg] ⁺	23.984493	1	pos
4	[M+K] ⁺	38.963158	1	pos
5	[M+Fe] ⁺	55.934388	1	pos
6	[M+Cu] ⁺	62.929049	1	pos
7	[M+2H] ⁺	2.015101	1	pos
8	[M+3H] ⁺	3.022926	1	pos
9	[M-H] ⁻	-1.007276	1	neg
10	[M+35Cl] ⁻	34.969401	1	neg
11	[M+Formate] ⁻	44.998203	1	neg
12	[M+Acetate] ⁻	59.013853	1	neg

The adducts library must have columns of `label`, `exact_mass`, `charge` and `ion_mode`.

We use this adducts file to adjust mass:

```
[8]: # if empty, use default adducts library
add_path = './data/adducts_short.tsv'
lib_add = anno.read_lib(add_path, ion_mode)
lib_add
```

```
[8]:
```

	label	exact_mass	charge
0	[M+H] ⁺	1.007276	1
1	[M+NH ₄] ⁺	18.033826	1
2	[M+Na] ⁺	22.989221	1
3	[M+Mg] ⁺	23.984493	1
4	[M+K] ⁺	38.963158	1
5	[M+Fe] ⁺	55.934388	1
6	[M+Cu] ⁺	62.929049	1
7	[M+2H] ⁺	2.015101	1
8	[M+3H] ⁺	3.022926	1

Now use function `comp_match_mass_add` to match compounds:

```
[9]: match_1 = anno.comp_match_mass_add(df, ppm, ref, lib_add)
      match_1
```

```
[9]:      id      mz molecular_formula      compound_name \
0      M152T40  152.043607      C5H8N2O2      Dihydrothymine
1      M154T37  154.062402      C8H8O3      p-Hydroxyphenylacetic acid
2      M154T37  154.062402      C8H8O3      3-Hydroxyphenylacetic acid
3      M154T37  154.062402      C8H8O3      ortho-Hydroxyphenylacetic acid
4      M154T37  154.062402      C8H8O3      Mandelic acid
5      M154T37  154.062402      C8H8O3      3-Cresotinic acid
6      M154T37  154.062402      C8H8O3      4-Hydroxy-3-methylbenzoic acid
7      M154T37  154.062402      C8H8O3      Vanillin
8      M157T35  157.036819      C4H10N2O2      2,4-Diaminobutyric acid
9      M157T35  157.036819      C4H10N2O2      L-2,4-diaminobutyric acid
10     M167T35  167.021095      C5H8N2O2      Dihydrothymine
11     M174T35  174.088395      C9H13NO      Phenylpropanolamine
12     M174T35  174.088395      C10H14O      Thymol
13     M174T35  174.088395      C10H14O      (S)-Carvone
14     M174T35  174.088395      C8H12O4      2-Octenedioic acid
15     M174T35  174.088395      C8H12O4      cis-4-Octenedioic acid
16     M181T36  181.060407      C8H8N2O3      Nicotinuric acid
17     M184T38  184.097942      C10H13N2      Nicotine imine
18     M185T39_2  185.082034      C5H15NO4P      Phosphorylcholine
19     M186T36  186.045606      C6H14N2O      N-Acetylputrescine
20     M187T38  187.097642      C5H15NO4P      Phosphorylcholine
21     M193T40  193.050761      C5H14N4      Agmatine
22     M200T36  200.061328      C7H16N2O      N-Acetylcadaverine
23     M201T39_1  201.051849      C10H10O3      4-Methoxycinnamic acid
24     M203T36_1  203.002108      C9H9NO      Indole-3-carbinol
25     M212T39  212.067866      C8H15NO3      Hexanoylglycine
26     M212T39  212.067866      C10H10O5      Vanilpyruvic acid
27     M217T37_1  217.018279      C10H11NO      Tryptophol
28     M221T37  221.012328      C9H11NO2      L-Phenylalanine
29     M223T38  223.008162      C4H10NO6P      O-Phosphothreonine
30     M223T40  223.096863      C12H14O4      Monoisobutyl phthalic acid
31     M226T44  226.128007      C8H18N4O2      Asymmetric dimethylarginine
32     M226T44  226.128007      C8H18N4O2      Symmetric dimethylarginine
33     M227T36  227.066175      C9H10N2O5      3-Nitrotyrosine
34     M229T38  229.069418      C4H10N3O5P      Phosphocreatine
35     M233T38  233.043479      C8H10N4O2      Caffeine
36     M245T44  245.045772      C7H15N3O3      Homocitrulline
37     M245T37_2  245.093315      C13H18O2      Ibuprofen
38     M249T38  249.038309      C8H10N4O3      1,3,7-Trimethyluric acid
39     M261T43  260.972975      C10H7NO4      Xanthurenic acid
40     M269T37_2  269.088048      C10H12N4O5      Inosine
41     M275T168  275.201932      C18H24O2      Estradiol
42     M275T168  275.201932      C18H24O2      17a-Estradiol
```

43	M277T181	277.217564	C18H28O2	19-Norandrosterone
44	M277T181	277.217564	C18H28O2	19-Noretiocholanolone
45	M278T71	278.148195	C11H20N2O6	Saccharopine
46	M279T233	279.233232	C18H30O2	alpha-Linolenic acid
47	M279T233	279.233232	C18H28O2	19-Norandrosterone
48	M279T233	279.233232	C18H28O2	19-Noretiocholanolone
49	M281T287	281.248903	C18H32O2	Linoleic acid
50	M281T287	281.248903	C18H30O2	alpha-Linolenic acid
51	M282T61	282.070271	C10H14N2O6	Ribothymidine
52	M282T61	282.070271	C10H14N2O6	3-Methyluridine
53	M283T37	283.103695	C11H14N4O5	1-Methylinosine

inchi \

```

0 InChI=1S/C5H8N2O2/c1-3-2-6-5(9)7-4(3)8/h3H,2H2...
1 InChI=1S/C8H8O3/c9-7-3-1-6(2-4-7)5-8(10)11/h1-...
2 InChI=1S/C8H8O3/c9-7-3-1-2-6(4-7)5-8(10)11/h1-...
3 InChI=1S/C8H8O3/c9-7-4-2-1-3-6(7)5-8(10)11/h1-...
4 InChI=1S/C8H8O3/c9-7(8(10)11)6-4-2-1-3-5-6/h1-...
5 InChI=1S/C8H8O3/c1-5-3-2-4-6(7(5)9)8(10)11/h2-...
6 InChI=1S/C8H8O3/c1-5-4-6(8(10)11)2-3-7(5)9/h2-...
7 InChI=1S/C8H8O3/c1-11-8-4-6(5-9)2-3-7(8)10/h2-...
8 InChI=1S/C4H10N2O2/c5-2-1-3(6)4(7)8/h3H,1-2,5-...
9 InChI=1S/C4H10N2O2/c5-2-1-3(6)4(7)8/h3H,1-2,5-...
10 InChI=1S/C5H8N2O2/c1-3-2-6-5(9)7-4(3)8/h3H,2H2...
11 InChI=1S/C9H13NO/c1-7(10)9(11)8-5-3-2-4-6-8/h2...
12 InChI=1S/C10H14O/c1-7(2)9-5-4-8(3)6-10(9)11/h4...
13 InChI=1S/C10H14O/c1-7(2)9-5-4-8(3)10(11)6-9/h4...
14 InChI=1S/C8H12O4/c9-7(10)5-3-1-2-4-6-8(11)12/h...
15 InChI=1S/C8H12O4/c9-7(10)5-3-1-2-4-6-8(11)12/h...
16 InChI=1S/C8H8N2O3/c11-7(12)5-10-8(13)6-2-1-3-9...
17 InChI=1S/C10H13N2/c1-12-7-3-5-10(12)9-4-2-6-11...
18 InChI=1S/C5H14N4O4P/c1-6(2,3)4-5-10-11(7,8)9/h4...
19 InChI=1S/C6H14N2O/c1-6(9)8-5-3-2-4-7/h2-5,7H2,...
20 InChI=1S/C5H14N4O4P/c1-6(2,3)4-5-10-11(7,8)9/h4...
21 InChI=1S/C5H14N4/c6-3-1-2-4-9-5(7)8/h1-4,6H2,(...
22 InChI=1S/C7H16N2O/c1-7(10)9-6-4-2-3-5-8/h2-6,8...
23 InChI=1S/C10H10O3/c1-13-9-5-2-8(3-6-9)4-7-10(1...
24 InChI=1S/C9H9NO/c11-6-7-5-10-9-4-2-1-3-8(7)9/h...
25 InChI=1S/C8H15NO3/c1-2-3-4-5-7(10)9-6-8(11)12/...
26 InChI=1S/C10H10O5/c1-15-9-5-6(2-3-7(9)11)4-8(1...
27 InChI=1S/C10H11NO/c12-6-5-8-7-11-10-4-2-1-3-9(...
28 InChI=1S/C9H11NO2/c10-8(9(11)12)6-7-4-2-1-3-5-...
29 InChI=1S/C4H10NO6P/c1-2(3(5)4(6)7)11-12(8,9)10...
30 InChI=1S/C12H14O4/c1-8(2)7-16-12(15)10-6-4-3-5...
31 InChI=1S/C8H18N4O2/c1-12(2)8(10)11-5-3-4-6(9)7...
32 InChI=1S/C8H18N4O2/c1-10-8(11-2)12-5-3-4-6(9)7...
33 InChI=1S/C9H10N2O5/c10-6(9(13)14)3-5-1-2-8(12)...
```

34 InChI=1S/C4H10N3O5P/c1-7(2-3(8)9)4(5)6-13(10,1...
 35 InChI=1S/C8H10N4O2/c1-10-4-9-6-5(10)7(13)12(3)...
 36 InChI=1S/C7H15N3O3/c8-5(6(11)12)3-1-2-4-10-7(9...
 37 InChI=1S/C13H18O2/c1-9(2)8-11-4-6-12(7-5-11)10...
 38 InChI=1S/C8H10N4O3/c1-10-4-5(9-7(10)14)11(2)8(...
 39 InChI=1S/C10H7NO4/c12-7-3-1-2-5-8(13)4-6(10(14...
 40 InChI=1S/C10H12N4O5/c15-1-4-6(16)7(17)10(19-4)...
 41 InChI=1S/C18H24O2/c1-18-9-8-14-13-5-3-12(19)10...
 42 InChI=1S/C18H24O2/c1-18-9-8-14-13-5-3-12(19)10...
 43 InChI=1S/C18H28O2/c1-18-9-8-14-13-5-3-12(19)10...
 44 InChI=1S/C18H28O2/c1-18-9-8-14-13-5-3-12(19)10...
 45 InChI=1S/C11H20N2O6/c12-7(10(16)17)3-1-2-6-13-...
 46 InChI=1S/C18H30O2/c1-2-3-4-5-6-7-8-9-10-11-12-...
 47 InChI=1S/C18H28O2/c1-18-9-8-14-13-5-3-12(19)10...
 48 InChI=1S/C18H28O2/c1-18-9-8-14-13-5-3-12(19)10...
 49 InChI=1S/C18H32O2/c1-2-3-4-5-6-7-8-9-10-11-12-...
 50 InChI=1S/C18H30O2/c1-2-3-4-5-6-7-8-9-10-11-12-...
 51 InChI=1S/C10H14N2O6/c1-4-2-12(10(17)11-8(4)16)...
 52 InChI=1S/C10H14N2O6/c1-11-6(14)2-3-12(10(11)17)...
 53 InChI=1S/C11H14N4O5/c1-14-3-13-9-6(10(14)19)12...

	inchi_key	exact_mass	adduct	ppm_error
0	NBAKTGXDIBVZOO-VKHYHEASA-N	152.04	[M+Mg] ⁺	3.52
1	XQXPVVBIMBYFF-UHFFFAOYSA-N	154.06	[M+2H] ⁺	-0.28
2	FVMDYYGIDFPZAX-UHFFFAOYSA-N	154.06	[M+2H] ⁺	-0.28
3	CCVYRRGZDBSHFU-UHFFFAOYSA-N	154.06	[M+2H] ⁺	-0.28
4	IWYDHOAUDWTVEP-ZETCQYMHSA-N	154.06	[M+2H] ⁺	-0.28
5	WHSXTWFYRGOBGO-UHFFFAOYSA-N	154.06	[M+2H] ⁺	-0.28
6	LTFHNKUKQYVHDX-UHFFFAOYSA-N	154.06	[M+2H] ⁺	-0.28
7	MWOOGOJBHIARFG-UHFFFAOYSA-N	154.06	[M+2H] ⁺	-0.28
8	OGNSCSPNOLGXSM-UHFFFAOYSA-N	157.04	[M+K] ⁺	-3.61
9	OGNSCSPNOLGXSM-VKHYHEASA-N	157.04	[M+K] ⁺	-3.61
10	NBAKTGXDIBVZOO-VKHYHEASA-N	167.02	[M+K] ⁺	-3.83
11	DLNKOYKMWOXYQA-VXNVDRBHSA-N	174.09	[M+Na] ⁺	-3.10
12	MGSRCZKZVOBKFT-UHFFFAOYSA-N	174.09	[M+Mg] ⁺	-3.23
13	ULDHMXUKGWMISQ-VIFPVBQESA-N	174.09	[M+Mg] ⁺	-3.23
14	BNTPVRGYUHFHJH-NWKNZROSA-N	174.09	[M+2H] ⁺	-1.52
15	LQVYKEXVMZXOAH-UPHRSURJSA-N	174.09	[M+2H] ⁺	-1.52
16	ZBSGKPYXQINNGF-UHFFFAOYSA-N	181.06	[M+H] ⁺	-1.99
17	GTQXYYYOJZZJHL-UHFFFAOYSA-N	184.10	[M+Na] ⁺	4.61
18	YHHSONZFQIEMCP-UHFFFAOYSA-O	185.08	[M+H] ⁺	4.80
19	KLZGKIDSEJWEDW-UHFFFAOYSA-N	186.05	[M+Fe] ⁺	3.25
20	YHHSONZFQIEMCP-UHFFFAOYSA-O	187.10	[M+3H] ⁺	4.52
21	QYPPJABKJHAVHS-UHFFFAOYSA-N	193.05	[M+Cu] ⁺	-0.69
22	RMOIHAKNOFHOE-UHFFFAOYSA-N	200.06	[M+Fe] ⁺	3.39
23	AFDXODALSZRGH-QPJXVBHSA-N	201.05	[M+Na] ⁺	-1.82
24	IVYPNXXAYMYVSP-UHFFFAOYSA-N	203.00	[M+Fe] ⁺	-3.42

25	UPCKIPHSXMXJOX-UHFFFAOYSA-N	212.07	[M+K] +	-2.29
26	YGQHQTMZRPHIBB-UHFFFAOYSA-N	212.07	[M+2H] +	-0.28
27	MBBOMCVGYCRMEA-UHFFFAOYSA-N	217.02	[M+Fe] +	-0.79
28	COLNVLDHVKWLRT-QMMMGPBSA-N	221.01	[M+Fe] +	-4.70
29	USRGIUJOYOXOQJ-GBXIJSLDSA-N	223.01	[M+Mg] +	-4.06
30	RZJSUWQGFCNFS-UHFFFAOYSA-N	223.10	[M+H] +	1.69
31	YDGMGEXADBOMJ-LURJTMIESA-N	226.13	[M+Mg] +	2.38
32	HVPFXCBJHIIJGS-LURJTMIESA-N	226.13	[M+Mg] +	2.38
33	FBTSQILOGYXGMD-LURJTMIESA-N	227.07	[M+H] +	-0.32
34	DRBBFCLWYRJSJZ-UHFFFAOYSA-N	229.07	[M+NH4] +	-0.94
35	RYYVLZVUVIJVGH-UHFFFAOYSA-N	233.04	[M+K] +	-0.23
36	XIGSAGMEBXLVJJ-YFKPBYRVSA-N	245.05	[M+Fe] +	0.17
37	HEFNWSXXWATRW-UHFFFAOYSA-N	245.09	[M+K] +	-2.13
38	BYXCFUMGEBZDDI-UHFFFAOYSA-N	249.04	[M+K] +	-0.56
39	FBZONXHGGPHHIY-UHFFFAOYSA-N	260.97	[M+Fe] +	4.13
40	UGQMRVRMYASKQ-KQYNXXCUSA-N	269.09	[M+H] +	0.01
41	VOXZDWNVPVJITMN-ZBRFXRBCSA-N	275.20	[M+3H] +	5.00
42	VOXZDWNVPVJITMN-SFFUCWETSA-N	275.20	[M+3H] +	5.00
43	UOUIARGWRPHDBX-CQZDKXCPSA-N	277.22	[M+H] +	4.90
44	UOUIARGWRPHDBX-DHMOVHTBWSA-N	277.22	[M+H] +	4.90
45	ZDGJAHTZVHVLLOT-YUMQZZPRSA-N	278.15	[M+2H] +	3.44
46	DTOSIQBPVRVQHS-PDBXOOCHSA-N	279.23	[M+H] +	4.93
47	UOUIARGWRPHDBX-CQZDKXCPSA-N	279.23	[M+3H] +	4.93
48	UOUIARGWRPHDBX-DHMOVHTBWSA-N	279.23	[M+3H] +	4.93
49	OYHQOLUKZRVURQ-HZJYTTRNSA-N	281.25	[M+H] +	4.97
50	DTOSIQBPVRVQHS-PDBXOOCHSA-N	281.25	[M+3H] +	4.97
51	DWRXFEITVBNRMK-JXOAFINSAN-N	282.07	[M+Mg] +	2.10
52	UTQUILVPBZEHTK-UHFFFAOYSA-N	282.07	[M+Mg] +	2.10
53	WJNGQIYEQLPJMNIOSLPCCCSA-N	283.10	[M+H] +	-0.01

1.4 Correlation Analysis

Next step is correlation analysis, based on intensity data matrix along all peaks. All results are filtered by the correlation coefficient, p-values and retention time difference. That is: keep correlation results in an retention time differences/window (such as 1 second) with correlation coefficient larger than a threshold (such as 0.5) and their correlation p-values less than a threshold (such as 0.05).

LAMP supports two correlation methods, `pearson` and `spearman`. Also parameter `positive` allows user to select only positive correlation results, otherwise positive and negative correlations will be used.

Two functions, `_tic` and `_toc`, record the correlation computation time in seconds.

```
[10]: thres_rt = 1.0
      thres_corr = 0.5
      thres_pval = 0.05
      method = "spearman" # "pearson"
```

```
positive = True
```

```
[11]: utils._tic()
corr = stats.comp_corr_rt(df, thres_rt, thres_corr, thres_pval, method,
                           positive)
utils._toc()
corr
```

Elapsed time: 0.7674360275268555 seconds.

```
[11]:
```

	name_a	name_b	r_value	p_value	rt_diff
0	M151T34	M153T34	0.80	1.267076e-23	0.02
1	M151T34	M155T34	0.71	1.752854e-16	0.20
2	M151T34	M161T34	0.78	1.869949e-21	0.14
3	M151T34	M163T34	0.69	3.239594e-15	0.20
4	M151T34	M167T35	0.51	5.776482e-08	0.73
...
1783	M283T34_1	M283T34_2	0.62	4.214876e-12	0.29
1784	M283T34_1	M285T34	0.82	5.937139e-26	0.08
1785	M283T34_2	M285T34	0.66	7.898957e-14	0.37
1786	M283T60	M284T60	0.86	1.033010e-29	0.15
1787	M283T339	M284T339	0.91	4.031333e-39	0.04

[1788 rows x 5 columns]

corr gives results of correlation coefficient(r_value), correlation p-values(p_value) and retention time difference(rt_diff).

Based on the correlation analysis, we can extract the groups and their sizes by:

```
[12]: # get correlation group and size
corr_df = stats.corr_grp_size(corr)
corr_df
```

```
[12]:
```

	name	cor_grp_size	cor_grp
0	M219T35	52	M221T34::M223T34::M225T35::M226T35::M229T34::M...
1	M215T35	52	M216T35::M217T35::M218T35::M219T34::M219T35::M...
2	M216T35	52	M217T35::M218T35::M219T34::M219T35::M221T34::M...
3	M217T35	52	M218T35::M219T34::M219T35::M221T34::M223T34::M...
4	M218T35	51	M219T34::M219T35::M221T34::M223T34::M225T35::M...
..
335	M196T593	1	M200T593
336	M272T59	1	M271T59
337	M257T51	1	M258T51
338	M223T54	1	M269T53
339	M271T59	1	M272T59

[340 rows x 3 columns]

1.5 Summarize Results

The final step gets the summary table in different format and save for the further analysis.

```
[13]: # get summary of metabolite annotation
sr, mr = anno.comp_summ(df, match)
```

This function combines peak table with compound matching results and returns two results in different formats. `sr` is single row results for each peak id in peak table `df`:

```
[14]: sr
```

```
[14]:
```

	name	mz	rt	exact_mass	ppm_error	\
0	M151T34	150.886715	34.152700	NaN	NaN	
1	M151T40	151.040235	39.838172	NaN	NaN	
2	M152T40	152.043607	40.303700	NaN	NaN	
3	M153T34	152.883824	34.174647	NaN	NaN	
4	M153T36	153.019474	35.785847	NaN	NaN	
..	
395	M283T61	283.068474	60.739869	NaN	NaN	
396	M284T108	284.223499	108.406389	NaN	NaN	
397	M284T339	284.267962	338.725056	NaN	NaN	
398	M284T60	284.195294	59.593561	NaN	NaN	
399	M285T34	284.775031	34.079641	NaN	NaN	

	molecular_formula	compound_name	inchi	inchi_key
0	NaN	NaN	NaN	NaN
1	NaN	NaN	NaN	NaN
2	NaN	NaN	NaN	NaN
3	NaN	NaN	NaN	NaN
4	NaN	NaN	NaN	NaN
..
395	NaN	NaN	NaN	NaN
396	NaN	NaN	NaN	NaN
397	NaN	NaN	NaN	NaN
398	NaN	NaN	NaN	NaN
399	NaN	NaN	NaN	NaN

[400 rows x 9 columns]

`mr` is multiple rows format if the match more than once from the reference file:

```
[15]: mr
```

```
[15]:
```

	name	mz	rt	molecular_formula	compound_name	inchi	\
0	M151T34	150.886715	34.152700	NaN	NaN	NaN	
1	M151T40	151.040235	39.838172	NaN	NaN	NaN	
2	M152T40	152.043607	40.303700	NaN	NaN	NaN	
3	M153T34	152.883824	34.174647	NaN	NaN	NaN	

4	M153T36	153.019474	35.785847		NaN		NaN	NaN
..	
404	M283T61	283.068474	60.739869		NaN		NaN	NaN
405	M284T108	284.223499	108.406389		NaN		NaN	NaN
406	M284T339	284.267962	338.725056		NaN		NaN	NaN
407	M284T60	284.195294	59.593561		NaN		NaN	NaN
408	M285T34	284.775031	34.079641		NaN		NaN	NaN

	inchi_key	exact_mass	ppm_error
0	NaN	NaN	NaN
1	NaN	NaN	NaN
2	NaN	NaN	NaN
3	NaN	NaN	NaN
4	NaN	NaN	NaN
..
404	NaN	NaN	NaN
405	NaN	NaN	NaN
406	NaN	NaN	NaN
407	NaN	NaN	NaN
408	NaN	NaN	NaN

[409 rows x 9 columns]

Now we merges single format results with correlation results:

```
[16]: # merge summery table with correlation analysis
res = anno.comp_summ_corr(sr, corr_df)
res
```

[16]:	name	mz	rt	exact_mass	ppm_error	\
0	M167T35	167.021095	34.882147	167.02	-4.57	
1	M276T36	276.077397	36.385373	276.08	-2.16	
2	M154T37	154.062402	37.183625	154.06	-3.84	
3	M181T36	181.060407	35.734801	181.06	2.39	
4	M174T35	174.088395	35.001130	174.09	-4.67	
..	
395	M279T50	279.159930	50.055451	NaN	NaN	
396	M279T79	279.163910	78.758079	NaN	NaN	
397	M282T85	282.207859	84.719202	NaN	NaN	
398	M283T47	283.110871	46.822069	NaN	NaN	
399	M284T108	284.223499	108.406389	NaN	NaN	

	molecular_formula	compound_name	\
0	C7H5NO4	Quinolinic acid	
1	C10H16N2O5S	Biotin sulfone	
2	C8H10O3	Hydroxytyrosol	
3	C6H7N5O2	8-Hydroxy-7-methylguanine	
4	C8H14O4	Suberic acid	

..
395	NaN	NaN
396	NaN	NaN
397	NaN	NaN
398	NaN	NaN
399	NaN	NaN

		inchi \
0	InChI=1S/C7H5NO4/c9-6(10)4-2-1-3-8-5(4)7(11)12...	
1	InChI=1S/C10H16N2O5S/c13-8(14)4-2-1-3-7-9-6(5-...	
2	InChI=1S/C8H10O3/c9-4-3-6-1-2-7(10)8(11)5-6/h1...	
3	InChI=1S/C6H7N5O2/c1-11-2-3(9-6(11)13)8-5(7)10...	
4	InChI=1S/C8H14O4/c9-7(10)5-3-1-2-4-6-8(11)12/h...	
..	...	
395		NaN
396		NaN
397		NaN
398		NaN
399		NaN

	inchi_key	cor_grp_size \
0	GJAWHXXKYYXBSV-UHFFFAOYSA-N	25.0
1	QPFQYMONYBAUCY-ZKWXMUHSA-N	13.0
2	JUUBCHWRXWPFFH-UHFFFAOYSA-N	12.0
3	VHPXSVXJBWZORQ-UHFFFAOYSA-N	9.0
4	TYFQFVWCELRYAO-UHFFFAOYSA-N	9.0
..
395	NaN	NaN
396	NaN	NaN
397	NaN	NaN
398	NaN	NaN
399	NaN	NaN

	cor_grp
0	M171T34::M197T36::M209T34::M211T34::M213T34::M...
1	M277T36_2::M278T36::M173T36_2::M186T36::M187T3...
2	M155T38::M158T37_2::M164T36::M171T37_2::M173T3...
3	M224T36::M225T35::M226T35::M227T36::M269T37_2::...
4	M211T34::M213T34::M219T34::M221T34::M229T35::M...
..	...
395	NaN
396	NaN
397	NaN
398	NaN
399	NaN

[400 rows x 11 columns]

The result data frame `res` is re-arranged as four parts from top to bottom:

- 1st part: identified metabolites, satisfied with correlation analysis
- 2nd part: identified metabolites, not satisfied with correlation
- 3rd part: no identified metabolites, satisfied with correlation
- 4th part: no identified metabolites, not satisfied with correlation

The users should focus on the first part and perform their further analysis.

You can save all results in different forms, such as text format TSV or CSV. You can also save all results into a `sqlite3` database and use [DB Browser for SQLite](#) to view:

```
[17]: f_save = False          # here we do NOT save results
      db_out = "test.db"
      sr_out = "test_s.tsv"

[18]: if f_save:
      # save all results into a sqlite3 database
      conn = sqlite3.connect(db_out)
      df[["name", "mz", "rt"]].to_sql("peaklist",
                                     conn,
                                     if_exists="replace",
                                     index=False)

      corr_df.to_sql("corr_grp", conn, if_exists="replace", index=False)
      corr.to_sql("corr_pval_rt", conn, if_exists="replace", index=False)
      match.to_sql("match", conn, if_exists="replace", index=False)
      mr.to_sql("anno_mr", conn, if_exists="replace", index=False)
      res.to_sql("anno_sr", conn, if_exists="replace", index=False)

      conn.commit()
      conn.close()

      # save final results
      res.to_csv(sr_out, sep="\t", index=False)
```

1.6 End User Usages

For end users, LAMP provides two computation options: command line interface(CLI) and graphical user interface (GUI).

To use GUI, you need to open a terminal and type in:

```
$ lamp gui
```

To use CLI, open a terminal and type in command with required arguments, something like:

```
$ lamp cli \
  --input-data "./data/df_pos_3.tsv" \
  --sep "tab" \
  --col-idx "1, 2, 3, 4" \
  --add-path "" \
```

```

--ref-path "" \
--ion-mode "pos" \
--cal-mass \
--thres-rt "1.0" \
--thres-corr "0.5" \
--thres-pval "0.05" \
--method "pearson" \
--positive \
--ppm "5.0" \
--save-db \
--save-mr \
--db-out "./res/test.db" \
--sr-out "./res/test_s.tsv" \
--mr-out "./res/test_m.tsv"

```

For the best practice, you can create a bash script `.sh` (Linux and MacOS) or Windows script `.bat` to contain these CLI arguments. Change parameters in these files each time when processing new data set.

For example, there are `lamp_cli.sh` and `lamp_cli.bat` in <https://github.com/wanchanglin/lamp/tree/master/examples>. You can run them and check the results in directory `examples/res`:

- For Linux and MacOS terminal:

```

$ chmod +x lamp_cli.sh
$ ./lamp_cli.sh

```

- For Windows terminal:

```

$ lamp_cli.bat

```